





ARTG ID : 427686

Australia Sponsor  
Australia Health Products  
Central Pty Ltd

# COVID-19/RSV/ Influenza A&B Antigen Test Kit

4-in-1 combination test kit,  
detect RSV, Influenza A/B and Covid-19 all at once.



## Methodology

Using the principle of double antibody sandwich method.



## Operation

Simple operation, easy to interpret,  
sampling and testing can be done  
anytime.



## Detection

The detection is fast, the result can be  
interpreted in 15 minutes.



## Accuracy

High sensitivity and specificity.



# Operation steps and result interpretation

**1**

Wash your hands.

**2**

Tear the aluminum foil on the extraction buffer tube. Place extraction tube into box tube stand.

**3**

Open the swab package and take out the swab.  
**Note: Do not touch the swab tip with finger.**

**4**

Tilt your head back slightly. Insert the swab about 1.5 to 2.5 cm into one nostril. Gently rotate the swab at least five times against the nasal wall.

**5**

Insert the same swab about 1.5 to 2.5 cm into the second nostril. Again, gently rotate the swab at least five times against the nasal wall.

**6**

Place the swab in the extraction tube and rotate the swab against the walls of the tube 5 times. Allow the swab to stand in the extraction buffer tube for 1 minute.

**7**

Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.

**8**

Press the nozzle cap tightly onto the tube.

**9**

Open the foil pouch and take out the test device.

**10**

3 drops must be added to the three specimen wells.

**11**

Read the result at 15 minutes. Do not interpret the result after 20 minutes.

**12**

Place the used test components into a biohazard specimen bag which can be sealed and dispose of in general waste.

**13**

Wash your hands thoroughly after test completion.

Instruction Video

Components	1 Test/Kit	5 Tests/Kit	20 Tests/Kit	25 Tests/Kit
1. Test Cassette	1 x	5 x	20 x	25 x
2. Extraction Buffer Tube	1 x	5 x	20 x	25 x
3. Disposable Swab	1 x	5 x	20 x	25 x
4. Biohazard Specimen Bag	1 x	5 x	20 x	25 x
5. Instruction for Use	1 x	1 x	20 x	25 x

Materials required but not provided : Timer

**20 independent 1 test/kit**

**COVID-19**

Positive Negative Invalid

**RSV**

Positive Negative Invalid

**Influenza A&B**

Positive Negative Invalid

**COVID-19 POSITIVE:** Two colored lines appear on the membrane. One line appears in the control region (C) and the other line appears in the test region (T).

**COVID-19 NEGATIVE:** Only one colored line appears in the control region (C). No colored line appears in the test region (T).

**COVID-19 INVALID:** Control line fails to appear.

**RSV POSITIVE:** Two colored lines appear on the membrane. One line appears in the control region (C) and the other line appears in the test region (T).

**RSV NEGATIVE:** Only one colored line appears in the control region (C). No colored line appears in the test region (T).

**RSV INVALID:** Control line fails to appear.

**Influenza A POSITIVE:** It is positive for Influenza A antigen if two Red lines appear. One red line should be in the control line region (C), and the other one appears in the A test line region.

**Influenza B POSITIVE:** It is positive for Influenza B antigen if two Red lines appear. One red line should be in the control line region (C), and the other one appears in the B test line region.

**Influenza A and B POSITIVE:** It is positive for both the antigens of Influenza A and Influenza B if three red lines appear. One Red line should be in the control line region (C), and another two should appear in A test line region and B test line region.

**NEGATIVE:** One Red line appears in the control region (C). No red line appears in the influenza A and B test region (T).

**INVALID:** Control line fails to appear.

## Product information

Product name	Specimen	Format	Pack/box	Shelf life	Storage temperature	Certificate
COVID-19/RSV/ Influenza A&B Antigen Test Kit	Nasal swab	Cassette	1T/5T/20T/25T	24months	2-30°C	TGA