

Helping you to protect what counts most

**3M™ Adaptic™ Touch Non-Adhering
Silicone Dressing**



Protects newly formed tissue

3M™ Adaptic™ Touch Non-Adhering Silicone Dressing is a flexible, open mesh silicone wound contact layer comprised of cellulose acetate coated with a soft tack silicone that has been designed to optimise fluid management and minimise adherence and pain at dressing change.

Protecting new formed granulation tissue is not only important to support successful wound healing, but also helps to protect the patient from pain and trauma which might be experienced at dressing change.

In a randomised control trial (RCT) on 70 patients, with burns, traumatic and post-surgical wounds 3M™ Adaptic™ Touch Non-Adhering Silicone Dressing scored as well as Mepitel® in terms of adherence, rates of wound closure, and optimal effectiveness.¹

Adaptic Touch Dressing supports the protection of the wound bed and its delicate wound tissue in three ways:



1

Providing wound bed protection and patient comfort

Adaptic Touch Dressing protects the wound bed through its gentle silicone contact layer and optimal pore size which are designed to minimise adherence to the wound bed and secondary dressing. It can be left in place for several days to protect vulnerable wound tissue and minimises patient pain at dressing change.²



2

Performance

The open mesh structure of Adaptic Touch Dressing allows free passage of exudate into an appropriate secondary dressing, reducing the risk of exudate pooling and maceration.³ The optimal pore size and transparent nature of the dressing allows visibility of the wound bed without the need to remove the dressing, providing undisturbed wound healing.³ The dressing is easy to use, does not stick to gloves, and remains in place *in situ* during dressing change.³



3

Providing atraumatic removal for patients

Its conformability and cuttability allows clinicians to dress difficult shaped wounds providing full protection for the wound bed, without the dressing unravelling or shedding fibres.⁴

When to use

3M™ Adaptic™ Touch Non-Adhering Silicone Dressing is indicated for use in the management of dry to highly exuding, partial and full-thickness wounds, such as:

- Venous leg ulcers
- Pressure injuries
- Diabetic foot ulcers
- Skin tears
- Surgical wounds
- Donor sites
- First and second degree burns



Adaptic Touch Dressing and Negative Pressure Wound Therapy (NPWT)

In a series of international case studies, Adaptic Touch Dressing was used as an alternative to traditional wound contact layers under NPWT to benefit patients across a range of wound types. Over the course of treatments, clinicians reported that Adaptic Touch Dressing was easy to use and remove and did not not adhere or cause pain during wear or upon removal.⁴ In addition, the dressing was shown to allow exudate to pass freely into the canister, without adhering or causing pain.⁴

Adaptic Touch Dressing at a glance:



**Protects
the wound
bed**



**Minimises
pain and
trauma**



**Can be
used in
conjunction
with NPWT**



**Suitable
for use on
fragile skin**



Cutttable

Ordering information

3M™ Adaptic™ Touch Non-Adhering Silicone Dressing

Code	Size	Dressings per box
TCH501	5 cm x 7.6 cm	10
TCH502	7.6 cm x 11 cm	10
TCH503	12.7 cm x 15 cm	10
TCH504	20 cm x 32 cm	5

To learn about the benefits of 3M™ Adaptic™ Touch Non-Adhering Silicone Dressing please contact your local 3M representative.

Explore the 3M™ Adaptic™ Dressing Range

3M™ Adaptic™
Non-Adhering Dressing



3M™ Adaptic™ Digit
Non-Adhering Dressing



References:

1. Pierrefeu-Lagrange et al. A randomised controlled study evaluating the clinical benefits of a cellulose acetate mesh coated with a soft tack silicone in the management of acute wounds. Wounds. 2018; 30(4). 2. Ivins N., Harding K. A series of case studies using a non-adhering silicone dressing. Poster presented at: EWMA; May 25-27, 2011; Brussels, BE. 3. Bianchi J. et al. Consensus guidance for the use of Adaptic Touch™ non-adherent dressing. Wounds UK. 2011; 7(3). 4. International case studies. Case series evaluation: Adaptic Touch™ in conjunction with negative pressure wound therapy. London: Wounds International, 2015.

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products. Please consult a healthcare provider and product instructions for use prior to application. This material is intended for healthcare professionals.



Medical Solutions Division
3M Australia Pty Ltd
Building A, 1 Rivett Road
North Ryde NSW 2113

Phone: +1300 363 878
Web: 3M.com.au/medical

**KCI Medical Australia
PTY Limited**
3M Australia Pty Ltd
Building A, 1 Rivett Road
North Ryde NSW 2113

Phone: +1300 524 822

InterMed Medical Limited
71 Apollo Drive
Albany, Auckland 0632
New Zealand

Phone: +09 415 4800

© 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorised use prohibited. ANZ-70-2013-1019-3 (12/21) AV011502943