clinell™

‘clean between’

to break the chain of infection



Have you missed out on Clinell environmental decontamination training? Catch up now.

Scan the QR code to learn how to protect patients from infections with Clinell Universal Wipes.





Wipe in an ‘S’ shaped pattern

Work from top to bottom

Wipe from clean to dirty

Ensure correct contact time

One wipe, one surface

GAMA Healthcare Australia Pty Ltd.,
Suite 1, 33-37 Duerdin St, Notting Hill VIC 3168, Australia.
+61 (0)3 9769 6600 | info@gamahealthcare.com.au | www.gamahealthcare.com.au

GHA230232

gama healthcare