Universal Wipes
Choosing the right product

clinell™ SURFACE

Detergent or disinfectant?

Detergent wipes

Best for cleaning surfaces by removing:
- Dirt and grime
- Dust
- Organic matter

Use when:
Routine cleaning is all that's required.
Use before disinfectant only wipe.



(Pack of 215 - CDW215AUS)

Disinfectant wipes

Best for disinfecting surfaces by reducing harmful microorganisms. TGA approved Class IIB disinfectant for cleaning and disinfecting shared medical devices.

Use when:
- Cleaning high-touch surfaces, shared equipment and examination spaces
- When contamination is suspected



(Pack of 200 - CW200AUS)

What is contact time – and why does it matter?

Contact time is the time a surface needs to remain wet for the disinfectant to work effectively. Wiping a surface dry before the stated contact time means the product may not have enough time to achieve its intended efficacy.

Always follow the manufacturer's instructions for use, including contact time and surface suitability.

Why choose a 2-in-1 wipe?

A combined detergent and disinfectant product cleans and disinfects a surface in a single step. For a busy practice, this has practical benefits:

- Cleans and disinfects in one step
- Simplifies workflows
- Reduces the need for multiple products
- Supports consistent cleaning practices
- Improves compliance

GAMA Healthcare Austrailia Pty Ltd.,
Suite 1, 33-37 Duerdin St,Notting Hill VIC 3168, Austrailia.
+61 (0)3 9769 6600 | info@gamahealthcare.com.au | www.gamahealthcare.com.au

GHA260248

gama healthcare