

cobas[®] b 101 system

Improve management of diabetes and dyslipidaemia with the HbA1c & Lipid monitoring PoCT system

The cobas b 101 system

The cobas b 101 includes tests for glycated haemoglobin (HbA1c) and lipids on the same device, thereby addressing a diagnostic need at Point of Care.

Main benefits

Ease of use and convenience:

- On the spot results in just 3 simple steps
- Adherent to all relevant standards and methods – ensuring test precision and full compliance to guideline requirements
- Direct sample application onto a disc from a finger prick of capillary or venous whole blood
- No need for tubes or pipettes for sample collection or transfer
- Individual or dual testing option
- Fast turnaround time – an intuitive 15 minute workflow from patient preparation to result of both HbA1c and Lipid Panel
- Increased convenience as both colour coded test discs can be stored at room temperature

Workflow efficiency:

- Patient data storage with 5,000 patient results, 500 QC results, and 50 Operator IDs
- Service free and no calibration is needed; samples and discs are checked for integrity, and all steps of the process are controlled
- QC menu features configurable test intervals and target ranges for reliable quality results
- IT connectivity solutions with cobas IT 1000 and printer option

\$3,875

RDS6378668

[Ask for your demo](#)

Reagents

cobas b 101 HbA1c

RDS6378676 Box/10 \$128.70

cobas b 101 Lipid

RDS6380115 Box/10 \$128.70



Compliance and effective monitoring...

...two essential components for optimising oral anticoagulation therapy

CoaguChek® XS Plus professional system

The professional device that provides the TOTAL SOLUTION for your practice

WHY COMPROMISE?
TOP PRODUCT

- Only 8µL drop of capillary blood from a finger prick
- Results in one minute enables immediate dosage adjustment
- "Out of meter" dosing ensures ease of use and avoids cross-contamination
- Quality Control: Each strip is checked to ensure integrity of each test
- Touch screen operations and icons
- Holds Patient ID, Operator ID & QC ID for quality compliance
- Large memory capacity: 1000 patient and 500 QC results can be viewed on screen
- Insensitive to heparin means rapid transition from heparin to warfarin therapy
- Unaffected by therapeutic levels of UFH/LMWH
- Patient record history: retrieve results by list
- Battery and/or mains operated
- Test can be done while holding unit in hand. No need to be placed flat
- Connectable to IT systems
- Includes Box/24 Strips
- **Lifetime Warranty**

CoaguChek®
Because it's my life



RDS58646401

\$850

Roche is the leading provider of INR devices in Australia and throughout the world*



INR online is a comprehensive anticoagulant management system (INR online Ltd) to help you monitor your patients on warfarin. It is available on the internet and can be used with both CoaguChek® XS systems in your practice.[†]

*Data on file

[†]Roche Diagnostics Australia takes no responsibility for the performance of INR online

Coaguchek® XS Practice/patient system

Ideal for selected patients testing at home or on home visits

All features as above but excludes:

- No Touch screen operation
- Memory capacity; Patient ID
- IT connection
- Mains operation (battery only)

RDS48115001

\$695



Coaguchek® Strips



RDS46253580	Box/24	\$169.20*
RDS46253150	Box/48	\$249*
	3+	\$232*
	6+	\$219*



Accutrend® Plus

Four tests in one instrument!

- Fast, safe, accurate - results in less than 3 minutes!
- No fuss – requires only single drop of capillary or venous blood
- “Out of meter” strip dosing and no cross contamination
- Tests Cholesterol, Glucose, Triglycerides and Lactate

RDS96418001 **\$285***

Accutrend® Plus Strips

Cholesterol Strips	RDS18262171	Box/25	\$95.43*
Glucose Strips	RDS47475173	Box/25	\$27.30*
Triglycerides Strips	RDS38144016	Box/25	\$87.12
BM-Lactate Strips	RDS12654016	Box/25	\$58.60



Accutrend® Control Solution

Cholesterol	RDS18289190	\$24.10*
Glucose	RDS84878190	\$17.12*
Triglycerides	RDS38152190	\$24.23
BM-Lactate	RDS47335190	\$49.52

ACCU-CHEK® Performa

ACCU-CHEK Performa

Blood glucose monitor



RDS61582014

\$49*

- Fast 5 second test time for quick and easy blood glucose testing
- Tiny sample size <0.6 µL with under-dosing detection. Up to 5 seconds to re-dose the same strip
- Fast fill strip. Blood sample is quickly absorbed by the test strip
- Compensates for temperature variables. Performs at 6°C to 44°C
- Compensates for the effects of humidity
- Expired test strip alert notifies user when strips are close to expiration date
- Ideal for medical centres



ACCU-CHEK Performa

Strips	RDS48874020	Box/100	\$52.26*
Control Solution	RDS61736001	Each	\$18.98*



TRUE Balance

Blood glucose monitor

- No coding
- Glucose oxidise strip technology
- 365 memory
- 14 and 30 day averages
- 10 second readings
- 1 microlitre of blood
- Software download
- Life time warranty

\$39.95*

NIP1959723



True Balance Strips

NIP1959731

Box/100 **\$59.95***



True Track Strips

(previous model)

NIP1948215

Box/50 **\$29.19***

Accu-Chek Safe-T-Pro Plus

There is a lot of care in this little lancing device. Specially designed for use in the healthcare environment.

- Single use lancet device
- 3 adjustable depth settings for comfort, convenience and hygienic assurance

Box/200 **\$68.67***
 2+ **\$64.24*** RDS03539200



Top Seller

Haemolance Plus®

Safety lancet for capillary blood sampling

Safe to use

Haemolance Plus is a single use product. The needle is hidden in the cartridge and retracts permanently after use.

Easy to use

The ergonomic shape and grip of the lancet makes it easy to use while giving confidence of a controlled, uniform and accurate sampling result.

Box/150 **\$73.40*** 2+ **\$69.00***



Low Flow
 • 25G needle size
 • 1.4mm puncture depth
 TER8HP7575



Normal Flow
 • 21G needle size
 • 1.8mm puncture depth
 TER8HP7568



Solofix Safety Lancets

Solofix® Safety single use lancet has been designed for modern capillary blood sampling

- Pre-activated safety mechanism
- Automatic irreversible needle retraction

Solofix Safety Fine

Low blood flow –for blood glucose monitoring

- Size 25G
- Puncture depth 1.5mm
- Blood sample < 10 µl



Box/200 **\$47.93*** BN6183010



Accu-Chek Active Strips	RDS44418186	Box/100	\$53.04*
Accu-Chek Active control solution	RDS46324195	Each	\$18.98*
Accu-Chek Sensor Comfort (Advantage) Strips	RDS49273621	Box/100	\$52.26*
Accu-Chek Sensor Comfort control solution	RDS30683001	Each	\$15.53*
Accu-Chek Go Strips	RDS82913386	Box/50	\$26.80*
CareSens Blood Glucose Test Strips	XXC008	Box/50	\$33.41*
Freestyle Optium Blood Glucose Test Strips	ADC9988470	Box/100	\$52.26*
Freestyle Lite Blood Glucose Test Strips	ADC7098470	Box/100	\$52.26*
Accu-Chek Softclix Lancets	RDS07484200	Box/200	\$28.01*
Accu-Chek Softclix Lancet Device	RDS07450001	Each	\$21.29*
Unistik 3 Normal Safety Lancets	CHAT1042	Box/100	\$42.24*
Unistik 2 Extra Safety Lancets	POIAT0712	Box/100	\$26.95*
Unistik 3 Extra Safety Lancets	POIAT1012	Box/100	\$34.95*
Ames Minilet Lancets	BAY5966	Box/200	\$74.49*

**More Strips,
 Control Solutions,
 Lancets &
 Lancet devices**

Can't find the lancet you use?
 Call Customer Service and we'll get it
 in stock for you.

