



\$249

MNL10122306

LUXO LFM 101 Magnifier

The 22 watt energy saving circular fluorescent lamp provides a high light output with a minimum of heat. The Ø 127mm glass lens (1.75 X magnification) is positioned in the middle of the circular tube; the ideal way of illuminating an object, virtually eliminating shadows. A hinged lens cover is standard; it acts like a dust cover and prevents the lens from acting like a burning glass.

- 5 year warranty
- Tube & clamp not included

22W Fluorescent Tube	RBKML235	\$35.00
AH Table Clamp (as in image)	MNL89900101	\$49.00



MAGGYLAMP

\$375

Wall Model **RBKML150WM**
Mobile Model **RBKML150MM**

Maggy-Vue Magnifying Lamp

- Top quality German design
- Circular lens, 120mm diameter, 2 X magnification
- 22 watt fluorescent lamp (included)
- Shadow-free diffused lighting
- Can be bench, wall or floor mounted
- Supplied with heavy duty bench clamp
- Spring balanced joints, easy to adjust
- Optional clip-on lens to create 3 X magnification



Maggy-Vue Clip on Lens	RBKML151	\$139.00
22W Fluorescent Tube	RBKML235	\$35.00
Mobile Base	REA001	\$149.00

Tevisio - LED magnifier light

With the latest LED's, innovative arm technology and a visual field designed for optimum performance, the TEVISIO magnifier luminaire combines high efficiency with ergonomic design.

Features:

- Unique lower arm assembly allows head to swing below bench/shelf height
- 48 premium SMD LEDs
- Continuous dimming to 10%
- 2 x magnification
- 6,000 lux output using only 14 watts
- Visualiser function for 3D shadow effects
- Large arm radius with innovative flexible head joint
- Balanced arm, easily adjustable
- Maintenance-free design with LED life up to 50,000+ hours
- Wide-angle lens Ø 153mm for distortion-free vision
- 784mm reach
- 2 year warranty

\$780

RBML160

The new LED MAGNIFIER from Waldman



Innovative flexible head joint

Continuous dimming to 10%

Can swing below bench or bracket height

Watch the Tevisio product video



SUPERLUX Task lights

- Height adjustable from 450mm to 580mm
- Even, diffused illumination
- Accepts incandescent or fluoro energy saver globes
- Globe not included
- Maximum reach 1070mm
- 1 year warranty

Desk/Mobile Model

- Unrestricted spring tension mechanism
- Ideal for mounting on mobile base

\$125

DAL1490

Wall Model

- Incorporated stop mechanism to prevent lamp damaging wall
- For wall mounting only

\$159

DAL1491



Desk Bracket

DAL1497

\$49.00

Mobile Base

REA001

\$149.00

Wall Bracket

DAL1496

\$28.00

Setting up a new surgery?

More and more doctors are choosing Team to set up their new practice.

We have the experience, the know how and of course the right price.

We'll guide you from start to finish.

Call us today!



ML-308

Halogen Exam light

The economical solution for general surgery lighting

- 70cm gooseneck
- Sure grip handle with on/off switch
- Warm white 2300°K colour
- 35 watt halogen 24° beam
- Wall or mobile mounted
- 4,000 hour bulb life
- 2 year warranty

\$295

RBKML308



MAGGYLAMP

Maggylamp Exam Accessories

Wall bracket RBKML307 \$65.00

Bi-fold wall bracket RBKML305 \$215.00

Mobile base bracket RBKML309 \$49.00

Mobile base REA001 \$149.00

ML-308 LED

LED Exam light

- 40,000 hour bulb life
- Cool white light, 5500°K colour
- Concise beam angle of 30°
- Low energy consumption, only 5 watt
- 3 x 1 watt LED giving 16,000 lux at 30cm
- Low heat output
- 2 year warranty



RBKML308LED

\$495

Bi-fold bracket

2 x 40cm arms giving 80cm reach from wall

Ideal for positioning lamp the full length of exam couch



Heine HL 5000® Examination Light

- Long-lasting 12V/50W halogen bulb gives 140,000 Lux light output at 30cm
- Extremely compact illumination-head. No obstruction of view, ease of handling
- 3 different mounting options: wall, table or wheeled stand
- Halogen bulb lasts at least 4000 hours
- Solid construction with high-specification 120cm radius gooseneck design
- Additional UV filter
- 5 year warranty



\$929

HEJ00527102

Examination Light HL 5000	Code	R.R.P	TEAM Price
HL 5000 exam light on lightweight wheeled stand L, polycarbonate base [05]	HEJ00527108	\$1288.42	\$1179
HL 5000 exam light with wall mount [04]	HEJ00527102 •	\$1006.15	\$929
HL 5000 exam light with clamp mount [03]	HEJ00527104	\$1166.28	\$1079
HL 5000 exam light spare bulb 12V/50W halogen [01]	HEJ00527111	\$63.65	\$59

• = popular

Heine HL 1200® Examination Light

- Long-lasting 12V/20W halogen bulb gives 90,000 Lux light output at 30cm
- Ultra-slim head for more freedom of movement and an unobstructed view
- 3 different mounting options: wall, table or wheeled stand
- Halogen bulb lasts at least 4000 hours
- Swivels though complete 180 degrees - flexible and easy to position
- Built in temperature filter avoids patient discomfort
- "Stay-put" construction of flexible section - adjust the light to virtually any position and it "stays-put," will not sag or droop
- Additional UV filter
- 5 year warranty



\$749

HEJ00527054

Examination Light HL 1200	Black	White	R.R.P	TEAM Price
HL 1200 exam light on lightweight wheeled stand L, polycarbonate base [05]	HEJ00527078	HEJ00527077	\$1100.42	\$1009
HL 1200 exam light with wall mount [04]	HEJ00527052	HEJ00527054 •	\$818.15	\$749
HL 1200 exam light with clamp mount [03]	HEJ00527056	HEJ00527058	\$978.28	\$899
HL 1200 exam light spare bulb 12V/20W halogen [02]	HEJ00527075		\$70.32	\$67

• = popular



[01]



[02]



[03]



[04]



[05]

Ordisi LED Examination light

The FLH-2 range is ideal for general and gynaecological examinations. Great quality examination lighting at an affordable price.

- 135,000 lux @30cm; 42,000 lux @50cm
- Lifetime LED, 50,000 hours
- Light field of 13cm Ø
- Colour temperature: 4,750°K
- Flexible gooseneck for optimum positioning
- Wall mount or mobile base
- 1 year warranty

Wall **ORDLH2131W**

\$525

Mobile **ORDLH2131M**

\$625

Call for your demo



Halux LED light

This LED examination lamp made by Derungs offers the best in LED medical lighting technology, combining low heat output with adjustable colour and brightness.

This gives the medical staff complete control when carrying out minor procedures and general examinations.

- Reach of 1 metre
- 5 light dimming levels
- Adjustable colour temperature
- Bright LED technology (no bulb change)
- 50,000 lux @0.5m
- Ceiling model with 1.2m bracket
- 2 year warranty

Mobile Model **RBKML505 + REA001** **\$1595.00**

Wall Model **RBKML503** **\$1475.00**

Ceiling Model **RBKML504** **\$1995.00**



Adjustable colour temperature – a big plus

The light colour temperature of these lamps can be changed by a simple touch of the control panel on the lamp head. The user can select the optimum colour for the examination depending on the skin area and type of tissue or wound. The three colour temperatures are selectable on the touch panel.

Adjustable colour temperature



warm white



neutral white



cold white

Also Available Halux Halogen Light

Mobile Model **RBKML500 + REA001** **\$1025.00**

Wall Model **RBKML501** **\$895.00**

Ceiling Model **RBKML502** **\$1349.00**

- 50 watt halogen beam
- Concise 10° beam and flexible head
- Ideal for examination/minor surgery
- 48,000 lux @0.5m
- 2 year warranty

UDP Examination Light

The UDP is an economical high intensity examination light widely used in consulting rooms and clinics. Ideal for general and gynaecological examination.


Strong and effective:

- Powerful even light beam of 11,200 lux at one metre
- Large task area action zone: 1900mm diameter x 1310mm height adjustment
- Cool beam enables extended procedures to be carried out without discomfort
- Lamp head has an aluminium inner reflector – self cooling effect for longer use
- All metal construction with powder coat finish for durability and ease of cleaning

Compact and flexible design:

- Arm mounted from above (2 metres from floor)
- Minimal workspace obstruction
- When not in use can be pushed up out of the way
- Time proven 'floating balance' suspension arm with adjustable fully enclosed springs
- Readily available lamp with average life of 4,000 hours; lamp change requires no tools
- Australian designed and manufactured
- Prompt supply of spare parts, technical service and support
- 2 year warranty

UDP Wall (with wall mount bracket and integral transformer)	MNLUDPUWTR	\$995
UDP Ceiling - ceiling mounted to a 2700mm floor to ceiling fixing point, with integral transformer	MNLUDPUCTR	\$995




ULED Examination Light

ULED is a significant step beyond traditional halogen examination lights. Specifically designed for close examination and minor procedure lighting tasks where a very powerful yet infinitely adjustable cold light is required.

- 20,000 lux at 1 metre
- 70 Lumens per watt; energy efficient
- Detachable head that can be easily replaced for ease of maintenance
- IP On/Off Switch and dimmable examination light
- Large task area action zone: Vertical: 1850mm Radius: 1050mm
- Arm mounted from above (2 metres from floor). When not in use can be pushed up out of the way
- Australian designed and manufactured
- 5 year warranty on led and 15 year warranty the support arm

ULED-UW Wall (with wall mount bracket and integral transformer)	MNLULEDUW	\$1,395
ULED-UC Ceiling - ceiling mounted to a 2700mm floor to ceiling fixing point, with integral transformer	MNLULEDUC	\$1,395



A key ULED feature is a push button to detach the light head from the arm and so facilitate upgrade, exchange or maintenance without a technician.





FlexLED

One of the most economical and flexible LED examination lights ever made in Australia.

Features

- 16,000 lux at 1m
- Colour temperature: 4,000K
- Colour rendering (Ra): 95
- Light beam angle: 10°
- Lux at working distance (approx 500mm) 64,000
- Hygienically sealed arm and head
- IP64 On/Off Switch
- Unlimited range of movement
- Total wattage including driver 12W
- Australian designed and manufactured
- 5 year warranty on LED
- 15 year warranty on arm

Available Jan 2014

MNLLLEDFLEX **\$695**

Desk Mount Clamp	MNLLLEDFLEXDMC	\$65
Desk Mount Permanent Screws	MNLLLEDFLEXDMP	\$65
Wall Mount	MNLLLEDFLEXWM	\$65
Mobile Base	MNLLLEDFLEXTM	\$299



LHH - Series

The LHH Series has widespread application throughout the health care sector; in hospitals, health centres and consulting rooms. It is one of the most extensively used examination lights in Australia.

It offers 3 distinct performance options; 2 Halogen versions and one outstanding LED version. The LHH series is designed to last the demands of today's clinical environment while requiring virtually no maintenance.

LHH 10 Designed for general examination; lights that require a large working light field

- Illumination: 4,430 lux at 1m
- Colour temperature: 3,100°K
- Colour rendering (Ra): 98
- Light Beam Angle: 24°


LHH 10 Wall	MNL8506 + MNL1181L	\$719
LHH 10 Ceiling	MNL8506 + MNLMEDIC600	\$1074

LHH 20 Designed for maximum illumination in a limited area and is ideal in areas like obstetrics, gynecology and wound care

- Illumination: 20,700 lux at 1m
- Colour temperature: 3,100°K
- Colour rendering (Ra): 98
- Light Beam Angle: 10°

LHH 20 Wall	MNL8507 + MNL1181L	\$865
LHH 20 Ceiling	MNL8507 + MNLMEDIC600	\$1220

5 Year Warranty



LHH LED

The LHH LED offers benefits over traditional halogen lamps. It is virtually maintenance free. The powerful LED light source provides an impressive 15,000 Lux at 1 metre and has an expected life span of 40,000 hours.

LHH is designed especially for demanding viewing conditions that require exceptionally good light such as medical examination and minor procedures.

- Colour temperature: 3,700°K
- Colour rendering (Ra): 93
- Light Beam Angle: 10°



LHH LED Wall	MNL8508 + MNL1181L	\$935
LHH LED Ceiling	MNL8508 + MNLMEDIC600	\$1290




Ordisi LED surgical light L21-25

The L21-25 range is ideal as a minor surgery light. Great quality theatre lighting at an affordable price.

- 100,000 lux @50cm; 50,000 lux @1m
- Lifetime LED, 50,000 hours
- 21cm diameter head
- Light field of 12cm Ø
- Colour temperature: 4,750°K
- Removeable sterilisable handle
- Flexible arm for optimum positioning
- Wall mount, ceiling mount or rollerbase (2 brakes)
- 1 year warranty

Wall [ORDL2125W](#)
 Ceiling [ORDL2125C](#)
 Mobile [ORDL2125M](#)

\$2,890



Centurion Minor Surgery Light

A popular minor surgery light. A proven performer for dermatology and plastic surgery.

- Elliptical reflector (380mm diameter) for good shadow reduction
- 43,000 Lux @ 1m
- 15cm light pattern diameter, with homogenous light distribution
- Colour correct light 4100° Kelvin
- Halogen bulb
- Smooth and stable articulating arm with full body coverage
- Removable, sterisable handle
- Made in the USA

Ceiling [DAL21544](#)

\$4,795

Phantom

The Phantom is equipped with a perfectly balanced floating arm which glides effortlessly at your fingertips. Its high performance LEDs produce 60,000 lux (at one metre) while consuming a mere 27-watts of total energy consumption.

Features

- 60,000 Lux at 1m
- Colour temperature: 4,000K
- Colour rendering (Ra): 95
- Light beam angle: 8°
- All metal construction
- On/Off Switch, dimmable Exam light
- Fully sealed head, that can be easily removed
- Arm has continuous rotation at 3 points
- 50,000 hours LED life over 15 -25 years in normal use
- Removable sterilisable handle
- Head shaped for ease of cleaning and sterilisation
- 15kg including ceiling mount
- Total wattage consumption 27W
- Australian designed and manufactured
- 5 year warranty on LED
- 15 year warranty on arm

Available Jan 2014



Phantom Wall
MNPLPHANW **\$3,795**

Phantom Ceiling
MNPLPHANC **\$3,795**



The Phantom is suitable for:

- Minor procedures
- Day surgeries
- Emergency areas
- Pre and post operative areas
- Intensive care



IP
64



team medical
supplies

Setting up a new surgery or needing an upgrade?

Lights are an important element of your surgery set-up.

Team Medical have the experience, the know-how and of course the right price.

Give us a call to discuss your lighting requirements.

We'll guide you from start to finish.