


Ansell

Hand Size Guide

Instructions

Place your right hand against the red line as shown. Read your hand size at the widest point from the scale shown opposite to determine your glove size.



Be Ansell sure.

Ansell



Micro Touch® Powder

Latex examination glove

- Powdered
- Natural rubber latex
- Textured surface with high grip



Box/100 \$8.42 10+ \$7.49 30+ \$6.68

Small AN4520N
 Medium AN4524N
 Large AN4526N



Micro Touch® DermaClean®

Latex examination glove

- Powder free
- Natural rubber latex
- Textured surface with high grip



Box/100 \$8.67 10+ \$7.65 30+ \$6.89

X-Small AN4570
 Small AN4572
 Medium AN4574
 Large AN4576
 X-Large AN4578



Micro Touch® Affinity®

Non sterile, Latex-free examination glove

- Powder free
- Nitrile formulation -Thiurams and mercaptobenzothiazoles free
- Textured surface with high grip
- Green coloured for easy identification as synthetic glove



Box/100 \$13.62 10+ \$12.76 30+ \$12.19

X-Small AN3770
 Small AN3771
 Medium AN3772
 Large AN3773



Micro Touch® NitraTex®

Nitrile examination glove

- Powder free
- Nitrile formulation - Thiurams free
- Textured surface with high grip
- Blue coloured for easy identification as synthetic glove



Box/100 \$13.91 10+ \$12.95 30+ \$11.99

Small AN4470
 Medium AN4472
 Large AN4474
 X-Large AN4476



Micro Touch® DentaGlove®

High quality latex dental examination glove

- Powder free
- Natural rubber latex
- Textured surface with high grip
- Surface finish



Box/100 \$11.75 10+ \$10.85 30+ \$10.25

Small AN4654N
 Medium AN4655N
 Large AN4656N



Micro Touch® DentaGlove® White Nitrile

High quality Latex-free examination glove

- Powder free
- Nitrile film has excellent puncture and chemical resistance
- Textured fingertips for high grip in wet or dry conditions
- 245mm glove length

X-Small AN4713
 Small AN4714
 Medium AN4715
 Large AN4716
 X-Large AN4717

Box/100 \$13.26 10+ \$12.55 30+ \$12.14



Gammex®

Powdered Latex Surgical Gloves

Provides an optimal combination of comfort, plus ease-of-use. It offers superior comfort and a fit that's natural to the hand. Its textured surface facilitates a better instrument grip in both damp and dry environments.

- 6 AN4151N
- 6.5 AN4152N
- 7 AN4153N
- 7.5 AN4154N
- 8 AN4155N
- 8.5 AN4156N



Box/40 \$48.64 6+ \$44.16 12+ \$40.96



Nutex® DermaShield®

Powder Free Sterile Surgical Gloves

Offers strength, softness, comfort and grip. The specific grip characteristics of Nutex DermaShield enable you to positively grip instruments and the smooth upper surface allows for effective manipulation of delicate items such as suture material.

- 6 AN4451D
- 6.5 AN4452D
- 7 AN4453D
- 7.5 AN4454D
- 8 AN4455D
- 8.5 AN4456D



Box/40 \$82.96 6+ \$77.58 12+ \$73.41



Gammex® PF

Powder-Free Latex Surgical Gloves

Gammex® PF latex surgical gloves provide excellent donning and comfort through the use of a hydrophilic/hydrophobic polymer lining allowing for both wet and dry donning. The wide thumb ball area combined with Ansell's soft latex formulation ensures greater comfort and less hand fatigue.

- Powder-free latex glove offering strength, softness, comfort and grip
- Mild textured finish for secure handling of instruments in wet and dry conditions
- Thinner and shaped finger pads for improved finger sensitivity

- 6 AN4151PF
- 6.5 AN4152PF
- 7 AN4153PF
- 7.5 AN4154PF
- 8 AN4155PF
- 8.5 AN4156PF



Box/40 \$88.00 6+ \$77.00



Gammex® PF Micro-Thin™

Powder Free Latex with Polymer Lining

Thin powder free gloves with polymer lining for increased sensitivity and precision. Gammex PF Micro-Thin are ideal microsurgery gloves.

- Powder-free
- Textured finger and palm area

- 6 AN4340PF
- 6.5 AN4341PF
- 7 AN4342PF
- 7.5 AN4343PF
- 8 AN4344PF
- 8.5 AN4345PF



Box/25 \$57.00 4+ \$52.65

Gammex® Non-Latex Accelerator-Free Sensitive

The new feel of advanced allergy protection featuring Sensoprene formulation

The new Gammex Sensoprene® formulation offers a comprehensive solution, protecting against both latex and chemical allergies. This unique surgical glove features unprecedented comfort and sensitivity as well as advanced allergy protection for healthcare workers and patients.

- Protects against both latex allergies and chemical allergies; Powder-Free
- Simplifies the glove ordering process – one glove to satisfy multiple requirements
- Designed for the ultimate sensitivity and comfort
- Allows for superior sensitivity while maintaining glove strength
- Protects against additional causes of skin reactions with features such as casein free formulation and CPC-free polymer coating
- Has the soft cream colour you want in a surgical glove

- 6 AN11481
- 6.5 AN11482
- 7 AN11483
- 7.5 AN11484
- 8 AN11485
- 8.5 AN11486



Box/25 \$79.86 4+ \$74.15

Gammex® PF DermaPrene®

Powder-Free Neoprene Surgical Gloves

Both latex-free and chemical accelerator-free. Ideal for those suffering from a latex allergy or chemical sensitivity. The Neoprene formulation produces a soft, comfortable glove with enhanced fit and feel.

- Powder-free and chemical accelerator-free
- Mild textured finish
- Green colour for instant identification as a synthetic glove

- 6 AN4211N
- 6.5 AN4212N
- 7 AN4213N
- 7.5 AN4214N
- 8 AN4215N
- 8.5 AN4216N



Box/25 \$81.00 4+ \$78.14

Flexi Powdered

Lightly powdered latex examination gloves

- Non sterile
- 240mm cuff
- Smooth



Small **MFXMLSH2**
 Medium **MFXMLSH3**
 Large **MFXMLSH4**

Box/100 **\$6.90** 10+ **\$6.05**

Flexi Powder Free

Powder free latex examination gloves

- Non sterile
- 240mm cuff
- Smooth



X-Small **MFXPFLH4**
 Small **MFXPFLH5**
 Medium **MFXPFLH6**
 Large **MFXPFLH7**
 X-large **MFXPFLH8**

Box/100 **\$7.30** 10+ **\$6.80**

Platinum Blue Nitrile

Powder free blue nitrile exam gloves

- Non sterile
- Latex free
- 240mm cuff
- Textured



Small **MFXPPLAT2**
 Medium **MFXPPLAT3**
 Large **MFXPPLAT4**

Box/150 **\$10.45** 10+ **\$9.20**

Cleantouch Vinyl

Powder free vinyl exam gloves

- Non sterile
- 240mm cuff



Small **MFXPFTGH1**
 Medium **MFXPFTGH2**
 Large **MFXPFTGH3**

Box/100 **\$5.45** 10+ **\$4.95**

Satin Pro Black

Powder free black latex exam gloves

- Non sterile
- 300mm cuff
- Textured



Small **MFXPSP1**
 Medium **MFXPSP2**
 Large **MFXPSP3**

Box/100 **\$12.85** 10+ **\$11.95**



Latex Examination Glove

Powdered



Small **AXP10115111**
 Medium **AXP10115112**
 Large **AXP10115113**

Box/100 **\$6.93** 10+ **\$6.32**



Latex Examination Glove

Powder free



Small **AXP10115121**
 Medium **AXP10115122**
 Large **AXP10115123**

Box/100 **\$7.55** 10+ **\$6.34**



Nitrile Examination Glove

Powder free

- Latex free



Small **AXP10115421**
 Medium **AXP10115422**
 Large **AXP10115423**

Box/100 **\$8.31** 10+ **\$7.60**

Protex® White Nitrile Glove

Not all nitrile gloves are the same. Whilst healthcare professionals require protection from exposure to infectious substances and cross-contamination, latex sensitivity and allergic reactions can be a concern. Protex White Nitrile gloves are manufactured from 100% high quality nitrile polymer and do not contain latex or fillers. They are manufactured in a facility that only produces nitrile gloves so there is no risk of contamination from other glove types.

- Excellent flexibility, stretch and comfort
- Superior puncture and tear resistance
- Excellent tactile sensitivity and dexterity
- Textured finger tips to enhance wet and dry grip

X-Small **SN36361431**
 Small **SN36361432**
 Medium **SN36361433**
 Large **SN36361434**
 X-large **SN36361435**

Box/100 **\$8.95**
 10+ **\$8.26**
 20+ **\$7.62**



Generation Pink® Examination Gloves

Show your support for breast cancer awareness - for every carton purchased Medline will donate \$1 to the National Breast Cancer Foundation.

- Latex free, advanced vinyl technology
- Powder free
- Designed for sensitivity, softness and comfort
- A great way to help in the fight against breast cancer

Small	MDLPINK6074			
Medium	MDLPINK6075	Box/100		\$6.16
Large	MDLPINK6076	10+		\$5.72



SensiCare Ice® Examination Gloves

- Non sterile Blue nitrile gloves
- Powder free, latex free
- Comfortable fit, reliable grip with excellent sensitivity

Small	MDL486801			
Medium	MDL486802	Box/200		\$16.20
Large	MDL486803	10+		\$14.40



Triflex® Surgical Gloves

- Sterile, latex, powdered
- Smooth anti slip finish
- Interlocking, beaded cuff to reduce roll-down

6	CAR2D7251	7.5	CAR2D7254		
6.5	CAR2D7252	8	CAR2D7255	Box/40	\$36.95
7	CAR2D7253	8.5	CAR2D7256	5+	\$34.97



Protexis® Surgical Gloves

- Latex with nitrile coating
- Powder free, sterile
- Durable with effortless donning
- Independent thumb design for flexible movement
- Interlocking, beaded cuff to reduce roll-down

6	CAR2D72N60X	7.5	CAR2D72N75X		
6.5	CAR2D72N65X	8	CAR2D72N80X	Box/50	\$70.07
7	CAR2D72N70X	8.5	CAR2D72N85X	4+	\$65.17



Protexis PI™ Surgical Gloves

With its unique synthetic polyisoprene formulation, Protexis PI provides a latex free alternative which offers the comfort, fit and feel of latex.

- Latex free - Synthetic polyisoprene
- Powder free, sterile
- Independent thumb design for flexible movement
- Interlocking, beaded cuff to reduce roll-down
- Replacement for Esteem gloves

6	CAR72PT60X	7.5	CAR72PT75X		
6.5	CAR72PT65X	8	CAR72PT80X	Box/50	\$120.06
7	CAR72PT70X	8.5	CAR72PT85X	4+	\$109.62



Choose latex-free for one less worry

Natural rubber latex allergy is a significant medical concern in healthcare today. Latex-sensitive patients and healthcare workers face a serious risk from any product containing latex, with exposure to latex gloves presenting a particular concern. To date, there is no known cure for latex allergy except eliminating exposure to latex products. Latex sensitivity can develop after repeated exposure to natural rubber latex.¹ The broken skin barrier caused by dermatitis increases the amount of exposure allowing latex proteins easy access through the skin.¹

When the level of sensitisation reaches the individual's unique threshold level, he or she will express symptoms on subsequent exposure to latex allergens and is considered latex allergic. Anyone with latex sensitivity is at risk of a life-threatening reaction and should be treated in the same manner as a latex allergic individual.²

Sensitisation develops in individuals genetically predisposed to latex allergy after multiple exposures to natural rubber latex over a highly variable period of time. The latency period ranges from several weeks to as long as 30 years.³

One of the quickest ways to reduce the risk of latex allergy is to switch to natural rubber latex-free gloves:

KIMBERLY-CLARK* STERLING* Nitrile Gloves


1. Sussman GL, Beezhold DH. Allergy to latex rubber. *Ann Intern Med* 1995;122:43-46.
2. AORN 2004 Standards, Recommended Practices, and Guidelines. "Latex Guideline" (2004):103, 106-108, 111-112.
3. Green-McKenzie A Hudes D. National Institute of Environmental Health Sciences. "Latex Induced Occupational Asthma in a Surgical Pathologist." *Environmental Health Perspectives* 113, no. 7 (July 2005):888.

Examination Glove Dispensers

- High-grade, clear plastic
- Can be mounted vertically or horizontally
- Fits most glove boxes



Single	KC50201	\$29.95
Double	KC50202	\$35.00
Triple	KC50203	\$41.50



Universal Glove Dispenser

- Powder coated metal
- Holds 3 boxes

Three tier
MIAGD101
\$49.00



Ansell Glove Dispenser

- Powder coated metal
- Fits Ansell gloves only

AN85 \$12.99





Synthetic Vinyl Glove

- Powder free
- Latex free - Polyvinyl Chloride
- 242mm cuff
- Smooth finish

Small KC55031
 Medium KC55032
 Large KC55033
 X-Large KC55034 (box/90)



Box/100 \$7.25 10+ \$6.50

Sterling Nitrile Glove

Fit, feel and tactile sensitivity of latex

- Powder free
- Latex free
- 242mm cuff
- Textured fingertip
- High box count

Small KC13490
 Medium KC13491
 Large KC13492



Box/200 \$19.25 10+ \$17.75

Try the KIMBERLY-CLARK* STERLING* Nitrile Exam Glove

We know you'll love them!

Why?

- A nitrile glove with superior barrier protection qualities
- Dexterity and tactile sensitivity of natural rubber latex
- Latex-free
 Unknowingly exposing a latex allergic patient to latex gloves and other latex products can be fatal.¹

How to choose the right size

1. Place the middle of your palm on the red dotted line on the sizing chart.
2. The height of your longest finger will be your correct glove size.

Request a sample

For a sample of KIMBERLY-CLARK* STERLING* Nitrile Exam Glove, email: Kimberly-ClarkHCAP@kcc.com

Offer valid until 31 April 2014.
 While stocks last.




1. Position Statement of the Organization of Main Nurse Executives. Approved January 26, 2001.

* Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. © 2013 KCWW.

L
LARGE

M
MED

S
SMALL



Sizing Chart



Kimberly-Clark

Trusted Clinical Solutions*