

Propulse NG Ear Irrigator



- Simple to operate and safer than traditional manual ear syringing procedures
- Operator and patient friendly
- Dual power supply; mains or rechargeable
- Quick release handle
- QRX system locks the jet tip into place during use
- Contrapulse water delivery system
- Variable, easily controlled water pressure
- Hands-free footswitch
- Soundproofed and double insulated motor for quieter operation
- Portable and rechargeable
- 12 month warranty

Items included:

- Propulse NG unit
- Hands-free footswitch
- 10 sterile disposable QRX tips
- 10 cleaning tablets
- Mains power adaptor
- Internal rechargeable battery
- User manual



\$419

BIOKIT4001

Propulse spares

Propulse QRX jet tips (for Propulse NG)	BIOKIT4002	Pkt/10	\$16.94
Propulse III jet tips (for Propulse III)	BIOKIT1001	Pkt/10	\$16.20
Propulse noots tank	BIOINS0021	Each	\$26.50
Propulse cleaning tablet	BIOCL0001	Box/200	\$35.00



Propulse® ProTect Cape

The Propulse® ProTect Waterproof Shoulder Cape helps protect the patient's clothes during the ear irrigation procedure.

The ProTect Cape can be used on either shoulder and is quickly and easily fitted to the patient. Its 3 ply design combines an absorbent paper material with a waterproof backing layer to provide maximum protection - the absorbent paper material minimises water "run-off" whilst the waterproof backing helps prevent water and earwax from seeping through to the patient's clothing.



BIOKIT1012
Box/36 \$17.50



Propulse® ProScoop Ear Wax Removal Scoop

The Propulse ProScoop has been designed to facilitate the manual removal of ear wax before, during and after the ear irrigation process.

- Scoop end designed by ear care specialist for easy removal of ear wax
- Indicator mark on scoop shaft at 15mm (±0.5mm) to show depth of penetration into ear
- Serrated end designed for easy application of cotton wool - for mopping the ear out after irrigation
- Non slip shaft for added safety



BIOMIR0028
Box/50 \$69



Ear Syringe

with two nozzles and shield, 100ml

BYDA2010 **\$46.50**



O ring for nozzle	BYDA2015	\$4.25
O ring for plunger	BYDA2012	\$4.95
Shield	BYDA2014	\$3.17
Nozzle Set	BYDA2013	\$19.45



NJ Philips Ear Syringe

Ear syringe NJ Phillips repair kit	MEAIKSPWX1187	Each	\$29.95
Ear syringe NJ Phillips nut	MEAITIS0001	Each	\$17.95
Ear syringe NJ Phillips nozzle only	MEAITIS020	Each	\$17.95
Ear syringe NJ Phillips silicone tips	MEASP1370B	Pkt/25	\$39.95

MEAITIS000
\$230



Suction Catheters Y-suction

6 FG	MMDAN042001
8 FG	MMDAN042002
10 FG	MMDAN042003
12 FG	MMDAN042004
14 FG	MMDAN042005
16 FG	MMDAN042006

Each **\$0.99**
50+ **\$0.85**



Suction Instrument Yankauer

Suction handle without vent
Sterile - double wrapped
CON30070025

Each **\$1.90**
50+ **\$1.75**

Suction handle with vent
Sterile - double wrapped
CON30900025

Each **\$1.95**
25+ **\$1.65**



Suction Tubing

2m non-sterile soft
MMD050000NS

Each **\$2.75**
10+ **\$2.25**

2m sterile double wrapped, soft
MMD050000

Each **\$1.70**
80+ **\$1.55**

3m sterile double wrapped, soft
MMD050001

Each **\$2.21**
25+ **\$2.02**



Disposable Suction Liners

The VacSax Advance range is a closed disposable suction system. The unique patented design of the bag, folded for ease of storage, will open immediately the vacuum is attached ensuring a perfect fit every time



1Ltr	ECO9910340	Each	\$5.27	25+	\$4.74
2Ltr	ECOVAL201	Each	\$4.49	30+	\$4.16

Manual Suction Pump

An effective and powerful manual suction pump for the quick suction of fluids from oro- and nasopharyngeal cavities of adults and children in emergency situations.

- Cost effective alternative to electrical suction pumps
- Lightweight and portable
- Easy to squeeze with maximum performance

ADTMS001PMP
\$149



Flaem Suction Pumps

Ideal for ENT, General Practice, Hospital Clinics, Aged Care & Homecare

- Ideal for immediate use
- Oil free, maintenance free
- Reliable double-piston mechanism
- Overflow protection (flotation shut off valve)
- Sturdy epoxy coated metal housing
- Antibacterial filter
- 12 month warranty

Flaem Suction Pro Pump
DALSP01P00

\$715



Flaem Port-a-Suction Pump
DALPS11P00

\$1025



Flaem Suction Pro Pump

- Maximum vacuum: -81 kpa/610mmHg
- Maximum flow: 20Lpm
- Jar size: 1000ml
- Overall dimensions: 330mL x 240D x 220H mm
- Weight: 5.5kg

Flaem Port-a-Suction Pump

- Fully portable battery powered unit which will operate with 240V power or 12 Volt DC
- Car adaptor included
- Maximum vacuum: -81kpa/610mmHg
- Maximum flow: 20Lpm
- Jar size: 1000ml
- Overall dimensions: 370L x 100D x 240H mm
- Weight: 5.6kg

CA-MI Askir 30 Suction Pump

The Askir 30 is a compact and easily portable surgical aspirator, ideal for use in a General Practice or at home. This suction pump is suitable for use with both adults and children.

- Designed for non continuous use – 120 min on / 60 min off
- Can be used for nasal, oral and tracheal aspiration of body fluids
- Comes with 1000ml autoclaveable jar with overflow valve system
- Dimensions: 350W x 210H x 180mm
- Suction pump 30 litre, 240 volt
- 2 year warranty



BIORE310100 **\$585**

DF-750 Portable Suction Unit

Ideal for General Practice

- Lightweight and easy to carry
- Low noise and powerful suction
- Overflow auto shut-off
- Twin diaphragm pump design
- 12 month warranty



- Mains operation
- 20 LPM Flow rate
- 650~700 mmHg max suction pressure
- 1L capacity
- Weight: 5.5kg
- Dimensions: 43L x 22W x 32H cm

MDQDF750 **\$599**

DF-506 Suction Unit

Ideal for General Practice, ENT, Day Surgery, Cosmetic Surgery

- Low noise
- Oil free, easy maintenance
- Overflow auto shut-off
- Manual or foot control operation
- Exceptional flow rate of 40 LPM
- Low power consumption
- 12 month warranty



- Mains operation
- Twin piston pump design
- 40 LPM Flow rate
- 720 mmHg max suction pressure
- 2L capacity
- Weight: 15.6kg
- Dimensions: 54L x 26W x 35H cm

MDQDF506 **\$1,559**

JX820D Portable Suction Unit

Lightweight and portable

- Provided with three kinds of electrical sources:
 - AC 100V-240V 50/60Hz
 - External DC 12V (car adaptor included)
 - Internal battery
- Can be conveniently used in many places
- Battery capacity indicator (green, yellow, red)
- Battery duration of full charge: ≥30min
- Long life piston pump design
- Overflow protection device
- 12 month warranty



- Mains or battery operated
- >20 LPM flow rate
- 600 mmHg max suction pressure
- 1L capacity
- Dimensions: 21L x 14W x 30H cm
- Weight: 4.5kg

MDQJX820D **\$625**