

Looking for a screening audiometer that saves you time?

The Entomed SA 202 performs automatic tests based on your input. It's a lightweight screening audiometer that will speed up your testing process with automatic tests based on your preferences saving you precious time.

SA202 – Screening Audiometer

A versatile stand-alone screening audiometer

SA 202 is a manual / automatic, true two-channel, pure-tone, air conduction screening audiometer with full frequency and output range. Its patented digital design and RoHS compatibility makes it well adapted to the modern, environmentally aware user. SA 202 is user friendly, robust, reliable and yet its light-weight design makes it ideal for frequent field use.

- Built-in talk forward microphone with volume control for convenient communication
- Internal memory - stores up to 50 complete audiograms in non-volatile memory
- Hard copy print-outs to any PCL 5 printer, via optional SA 07 printer interface cable
- Computer interface, via optional RS 232 computer cable

SA202 Audiometer

1 x EN100212 = Entomed SA202 screening audiometer
1 x EN200195 = DD45/AUDIOCUPS with 1/4" stereo jack

ECOEN100212KIT1

\$3495

SA202 Audiometer with Software

1 x EN100212 = Entomed SA202 screening audiometer
1 x EN200195 = DD45/AUDIOCUPS with 1/4" stereo jack
1 x EN101181 = Entomed AUDIMAX II+ Software

ECOEN100212KIT2

\$4695

** The SA 202 is a manual/automatic IEC 60645-1 Type 4 puretone, air conduction screening audiometer and meets or exceeds the following standards: IEC 60645-1; ISO 389-1; ISO 389-2; ISO 8253-1; IEC 60601-1 Medical Electrical Equipment Requirements for Safety and Performance IEC 60601-1-2 Medical Electrical Equipment Requirements for EMC Compatibility.



1YR
warranty



For more information download the product brochure from our website

AUDIOLYSER® ADL20 Electronic Audiometer

Simple and quick to use, the ADL20 Audiometer integrates the latest digital technology. Its DSP (digital sound processor) generates pure sounds while ensuring optimal precision. The ADL20 is a versatile, portable, computerised, new generation professional instrument. A vital instrument to accurately screen and check for auditory problems.

- Class 4 tonal audiometer
- Portable: compact and lightweight
- Multilingual Vocal Assistance
- Computerised: Audiowin®20 software
- Standards: EN 60601-1, EN 60601-1-2, EN 60645-1, ISO 8253-1, ISO 7029, ISO 389-1, ISO 389-5, ISO 389-8

Standard Accessories:


- 1x User Manual
- 1x Audiowin® 20 software
- 1x License per work station
- 1x Transport cover

ECOADL20

\$2895



2YR
warranty



For more information download the product brochure from our website

