

CPR ADVISOR



The HeartSine Samaritan 500P provides clear voice prompts, ensuring that the rescuer knows what to do at all times. The 500P will also monitor chest compressions and will prompt the rescuer when they need to push harder, faster or slower. This means that High Quality CPR will be guaranteed at all times.

During CPR, if rhythm slows due to fatigue, or is too fast, the Samaritan 500P will prompt **“Push Faster”** or **“Push Slower”**.

Chest compressions need to reach a depth equal to 1/3 of the patient’s chest (As per ANZCOR guidelines). If depth is not adequate, the 500P will prompt **“Push Harder”** ensuring enough oxygen reaches the brain and vital organs, dramatically increasing the chances of survival.

If the rhythm and depth of compressions are adequate, the 500P will reassure the rescuer with **“Good Compressions”**. This has been proven to be invaluable in a stressful situation to provide the rescuer with confidence and reassurance.

CPR Advisor is a patented technology unique to the HeartSine 500P. It is the **ONLY** AED to provide CPR feedback based on clinical evaluation of the patient.