



Make Every Minute Count

AFINION™ HbA1c

Consultations with patients who do not have their lab results are challenging and inefficient. Point-of-care HbA1c testing can improve operational efficiency and lead to increased patient motivation and satisfaction.¹⁻⁴



3-MINUTE TEST TIME



HIGHLY ACCURATE RESULTS



NO INTERFERENCE FROM COMMON Hb VARIANTS



NGSP AND IFCC CERTIFIED

AFINION™ HbA1c

The Afinion™ HbA1c is a simple fingerstick test that provides a patient's HbA1c value within minutes, during the consultation.

- 1.5 µL sample volume
- Capillary and venous whole blood sample material
- Measuring range: 4%–15% HbA1c/20–140 mmol/mol HbA1c
- 3 months of room temperature storage
- Ready-to-use controls in two levels

BENEFITS OF POINT-OF-CARE TESTING

The 3-minute Afinion HbA1c test fits perfectly within the workflow:

‘Point-of-care HbA1c testing represents an opportunity for improved care of people with diabetes.’

— Egbinike V, et al. *Diabetes Educator*. 2013;39:66-73.

‘Point-of-care testing for HbA1c provides the opportunity for more timely treatment changes.’

— ADA. *Diabetes Care*. 2018;41(Suppl 1):S55-S64.

IDF RECOMMENDATIONS

International Diabetes Federation (IDF) recommendations for the use of HbA1c in type 2 diabetes⁵:

- Diagnostic criteria: HbA1c \geq 6.5% or HbA1c \geq 48 mmol/mol
- Provide the HbA1c result, as measured at the site of care or in the laboratory, before the clinical consultation
- Measure HbA1c every 2 to 6 months depending on level, stability of blood glucose control and changes in therapy

PRODUCT NAME	PRODUCT CODE
AFINION™ HbA1c	1116795
AFINION™ ACR	1116781
AFINION™ 2 ANALYZER	1116770
AFINION™ CRP	1116787
AFINION™ LIPID PANEL	1116801

CONTACT US TODAY FOR MORE INFORMATION.

TO ARRANGE A DEMONSTRATION, VISIT [GLOBALPOINTOFCARE.ABBOTT](https://www.abbott.com/globalpointofcare)

T: 1800 622 642 (FREE CALL AU) | T: 0800 466 951 (FREE CALL NZ) | T: +61 7 3363 7100

RAPIDDX.ANZ.ENQUIRIES@ABBOTT.COM

1. Shepard MD. *Clin Biochem Rev*. 2006;27(3):161-70.
2. Laurence CO, Gialamas A, Bubner T. *Br J Gen Pract*. 2010;60(572):e98-e104.
3. Miller CD, Barnes CS, Phillips LS, et al. *Diabetes Care*. 2003;26(4):1158-63.
4. Crocker JB, Lee-Lewandrowski E, Lewandrowski N, et al. *Am J Clin Pathol*. 2014;142:640-6.
5. IDF. *Global Guideline for Type 2 Diabetes*. 2012. ISBN 2-930229-43-8.

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. A2910.3 12/20

3-STEP PROCEDURE

- 1 Collect the sample with the integrated sampling device.



- 2 Place the sampling device back in the test cartridge.



- 3 Place the test cartridge in the analyzer and close the lid. The processing starts automatically.

