

Dual Frequency Ultrasound Unit



Shown with optional
1cm Ultrasound Head



Nu-Tek® Ultra Rehab2 – UT2200

Dual Frequency Ultrasound with 5cm² Transducer

Features & Benefits:

- 1 & 3MHz ultrasound frequency
- Pulsed and continuous modes
- Variable duty cycles (16Hz, 48Hz and 100Hz)
- Ergonomically designed sound heads in 1cm², 5cm² sizes
- Contact control with both visual and audio feedback
- 20 pre-set programs & 80 custom programs provided

Screen Shot:



Specifications:

Power Supply:	AC100-240V, 47-63Hz
Intensity (Max):	3 W/cm ²
Ultrasound Frequency:	1 MHz/3MHz +/- 10%
Pulse Frequency:	16Hz, 48Hz, 100Hz
Duty:	10%, 20%, 30%, 40%, 50%, 60%, 70%, 80%, 90%
ERA:	1cm ² , 5cm ²
BNR:	< 5.0
Timer:	1 to 30 minutes
Dimensions:	285mm (L) x 197mm (W) x 153mm (H)
Weight:	2kg
Certification:	FDA clearance, EU MDD, CFDA, TGA, CMDCAS
Electrical Safety Class:	Class I, Type BF
Safety Test:	IEC 60601-1, IEC 60601-2-5, IEC 60601-2-10

Standard Accessories:

NUULTRAHEAD5	Nu-Tek Ultrasound Head 5cm
NUPOWERCORD	Nu-Tek Power Cable
ACGEL250	Ultrasound Gel - 250ml

Optional Accessories:

NUULTRAHEAD1	Nu-Tek Ultrasound Head 1cm
NUTROLLEY	Nu-Tek Trolley
NUBATTERY	Nu-Tek Battery Pack
ACGEL5L	Ultrasound Gel - 5 Litres

Code:

Description:

NUULTRASOUND	Nu-Tek® Ultra Rehab2 - UT2200
---------------------	-------------------------------

WHITELEY
HEALTHCARE

Phone: 1300 303 755 Fax: 1300 303 788

Email: sales@whiteleyhealthcare.com.au

Mail: PO Box 576, Hornsby NSW 1630 Australia



www.whiteleyhealthcare.com.au