



The AQ Medical 110 vaccine refrigerator is a compact, robust and cost effective unit, purpose-built for vaccine storage following the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines.

It was designed with low noise output and consistent temperature readings in mind, and maintains a stable temperature between the recommended +2°C to +8°C. The unit holds safety features such as visual and audible alarms in the event that temperatures fluctuate outside the set range.

All AQ Medical refrigerators are accredited by the Quality Care Pharmacy Program (QCPP), a quality assurance program that was developed by the Pharmacy Guild of Australia in consultation with the Pharmaceutical Society of Australia.

AQ Medical units are also backed by the extensive Avem Quirks service network. With qualified technicians and warehousing facilities across Australia, we can deliver responsive service and exceptional support on a national level.

The AQ Medical 110 vaccine refrigerator offers prime performance and quality at a competitive price.









## **AQ Medical 110 Vaccine Refrigerator**

## **Key Features**

- QCPP cold chain accredited
- Digital temperature control system preconfigured to hold temperature from +2°C to +8° C in line with national 'Strive for 5' storage guidelines
- · Audible and visual auto resetting alarms
- Adjustable shelves
- Lockable door
- · Low energy consumption
- Low noise output
- Dedicated technical service line

## **Optional Extras**

- Uninterruptible Power Supply unit with back up battery
- AQ Box cloud temperature monitoring system
- · Comprehensive programmed maintenance

## **Product Specifications**

Dimensions	Height 820 mm Width 595 mm
	Depth 595 mm
Gross Weight	41.5 kg
Gross Capacity	145 L
Temperature Range	+2°C to +8°C
Number of Shelves	4
Power Consumption/24 hrs.	1.26 kWh
Refrigerant	R600a
Warranty (parts and labour)	2 years

