


# AAXIS® Surgical Mask Type IIR

- High quality, 3-ply Surgical mask
- Complies with EN 14683: 2019 Type IIR
- ≥99.9% Bacterial Filtration Efficiency (BFE)
- Convenient dispenser box of 50
- Non-sterile

**Code:** **AXP10025519**

**Box/50**



## Test Data

Test Requirement	EN 14683:2019 Standard Type IIR	Aaxis Surgical Mask
Bacterial Filtration Efficient (BFE) %	≥98	Pass
Synthetic Blood Resistance	≥ 16.0 kPa >120 mmHg	Pass
Differential Pressure	<60 Pa/cm² <6.0 mmH2O/cm²	Pass
Microbial Cleanliness	≤30 cfu/g	Pass



**Looking after your team.** We are a proud independently owned Australian company

### Team Medical Supplies

Phone 1300 22 44 50 info@teammed.com.au  
Fax 1300 22 44 60 www.teammed.com.au

### BP Medical (WA)

Phone 1300 33 55 60 info@bpmedical.com.au  
Fax 1300 33 55 70 www.bpmedical.com.au



To view our terms and conditions please visit our websites - [www.teammed.com.au](http://www.teammed.com.au) or [www.bpmedical.com.au](http://www.bpmedical.com.au).

© Copyright Team Medical Supplies - August 2020. All Rights Reserved.