


SUREGARD Dual-Fitting Filter

Compatible with the Ndd Spirette and GANSHORN ScoutTube

The *Suregard Filter Mouthpiece* by Bird Healthcare is a low cost, infection control option for Ndd spirometry users. A new **direct fit Suregard filter** is available for the distal/exhaust end of the Ndd Spirette and Ganshorn ScoutTube, in addition to the mouthpiece end of the ScoutTube.



The *Suregard dual-fitting filter* for Ndd spirometers are **100% Australian manufactured**, utilising tec200 hydrophobic filter material which has been validated for filtration of bacteria and viruses at high (spirometry) flow rates. Filtration reports available on request.



GANSHORN ScoutTube

Compatible with the Ndd EasyOne and the Ganshorn SpiroScout

A recent **independent study** by Alexandra I. Smith, Jay R. Flack, Department of Respiratory and Sleep Medicine St George Hospital, Sydney, Australia, assesses the suitability of generic flow tube consumables and concluded that that Ganshorn ScoutTube do in fact **work with Ndd equipment, passing all ATS waveforms**. For access to the full independent study and for the manufacturers supporting evidence for the use of ScoutTube send a request to info@birdhealthcare.com.

Description	Code	TGA	QTY	Material Type	Recyclable
Suregard Mouthpiece filter	RJVKB14	332993	100	 Polypropylene	Yes
ScoutTube Mouthpiece	BHC-019420620-50	260378	50	 Polyethylene	Yes

Bundles Available

RJVKB14-ST-NC
100 Pink Pearl filters + 100 ScoutTubes + 100 Nose clips

Samples Available

RJVKB14-ST
100 Pink Pearl filters + 100 ScoutTubes

Bulk deals available

sales@birdhealthcare.com
www.birdhealthcare.com
1300 365 561



BIRD HEALTHCARE