



Proudly Supporting

Asthma
Australia

chamber **mini**

perfect for
laboratories
and training!

Actual Size

Butterfly Valve Technology™

offers low resistance to
exhaled air and directs the
flow away from the chamber

Anti-Static

eliminating the need
to rinse, drip dry or
prime prior to use

Compact & Low-Cost

low-cost spacer at a
compact size is ideal
for laboratory and
training purposes

Effective Delivery

provides effective drug
delivery for such a small
device, with ~30 mcg of
available medication



**BLISTER
PACKED**



Exclusive to
 BIRD HEALTHCARE

1300 365 561

www.birdhealthcare.com

SUPPLIED BY:

 