



The most trusted bacterial and viral filter  
on the Australian and New Zealand  
markets for over 2 decades



**BIRD HEALTHCARE**

[www.birdhealthcare.com](http://www.birdhealthcare.com)

**Key Features:**

- *Manufactured from inert, latex free, eco-friendly polypropylene plastic*
- *Utilises world leading filter material for superior filtration at high flow rates*
- *All in one integrated elliptical mouthpiece provides maximum seal, comfort and prevention of leaks*
- *Colour coded for easy identification between multiple systems*

**Material** *GVS (formerly Airsafety)*

**Filter Media** *TEC 200gm Electrostatic filter media. Electrostatic hydrophobic double outer scrim*

**Housing** *Polypropylene plastic (100% recyclable)*

**Mouthpiece Shape** *Oval, Elliptical*

**Mouthpiece Adapters** *Paediatric/small mouthpiece adapters available*

**Dead space** *50ml*

**Resistance** *0.71cmH<sub>2</sub>O/l/s at a flow of 12L/s*

**Filter efficiency** *Bacterial/Viral filtration 99.9999%*

**Compatible with** *Jaeger, Sensormedics, Medisoft, Ganshorn, MedGraphics, nSpire, Ndd, Vitalograph, Thorasys, CosMed, MIR.*

**Packaging** *100 per box (Grip nose clips optional)*

**Patented Ergonomic Mouthpiece**



**Paediatric Mouthpiece Adapter Also Available**



**Grip Nose Clip Also Available**



The SureGard® filter is an Australian owned, highly effective filter barrier device. The SureGard® provides both technicians and patients with optimal protection from cross contamination and infection. The R&D of the SureGard® is attributed to contributions from the medical and scientific community of Australia and New Zealand.



**RJVKB1 (Pink)**

Direct fit nSpire, MedGraphics (internal fit), MicroMedical, MIR SpiroBank, MIR MiniSpir, MIR Spiro Lab, MIR SpiroDoc Fukuda and Peak flow meters.  
External Diameter- 28.5mm  
Internal Diameter- 25mm



**RJVKB2 (Purple)**

Direct fit Viasys, SensorMedics, Jeager, QRS SpiroCard.  
External Diameter- Slight taper from 33mm to 32mm  
Internal Diameter- 29mm



**RJVKB3 (Blue)**

Direct fit MedGraphics (internal fit), MicroMedical, MIR Spiro Lab, MIR SpiroDoc, MIR SpiroBank, MIR MiniSpir, Peak flow meters, Viasys, SensorMedics, Jeager, QRS SpiroCard.  
External Diameter (right)- Slight Taper from 33mm to 32mm  
Internal Diameter (right)- 29mm  
External Diameter (left)- 28.5mm  
Internal Diameter (left)- 25mm



**RJVKB4 (Green)**

Direct fit Vitalograph Spirometers, MicroMedical, Spider USb, MicroLoop, SpiderXpert, MIR systems & CareFusion  
External Diameter- 32.5mm  
Internal Diameter- 29mm



**RJVKB5 (Red)**

Direct fit for all MediSoft systems  
External diameter- 31mm  
Internal diameter- 29mm



**RJVKB6 (Clear)**

Direct fit Ventilation systems, nebulisers, CPAP & BiPAP equipment.  
External Diameter (right)- 25.5mm  
Internal Diameter (right)- 22mm  
External Diameter (left)- 22mm  
Internal Diameter (left)- 19.5mm



**RJVKB7 (Pearl)**

Direct fit ndd EasyOne Spirometer (spirettes)  
External Diameter - 20mm  
Internal Diameter - 18.5mm



**RJVKB8 (Orange)**

Direct fit for all MedGraphics systems (external fit)  
External diameter- 36mm  
Internal diameter- 34mm



**RJVKB9 (White)**

Direct fit for all CosMed Pony, CosMed MicroQuark  
External diameter- 30mm  
Internal diameter- 28mm



**RJVKB10 (Pearl/Purple)**

Direct fit EasyOne Air (FlowTube)  
External width - 24.5mm



**RJVKB11 (Orange/White)**

Direct fit TremeFlo FOT  
External diameter- 36mm  
Internal diameter- 34mm



**RJVKB12 (White/Orange)**

Direct fit MedGraphics (Clip On)  
External diameter- 36mm  
Internal diameter- 34mm



[www.birdhealthcare.com](http://www.birdhealthcare.com)  
[sales@birdhealthcare.com](mailto:sales@birdhealthcare.com)

**Australia**  
Melbourne

**New Zealand**  
Auckland

**France**  
Paris

**England**  
London

**Hong Kong**  
Hong Kong