



# Portable Fetal Doppler



CE 0123

Model: A100D



## Feature

- Lightweight and Easy-To-Use
- 1.77Inch TFT Color LCD screen
- 2Mhz high sensitive probe
- Big font FHR display and FHR graphic display
- 500MAH rechargeable battery with direct charging function
- Real Time FHR and (3 minutes) FHR curve modes

## Specifications

<b>FHR</b>	Ultrasound Frequency:	2.0 MHz
	Intensity:	5mW/cm2 or less
	FHR Range:	60-210 bpm
	FHR Range for Manual Count mode:	30-250 bpm
	FHR Accuracy:	± 1 bpm

<b>Standard</b>	2.0 MHz Doppler probe:	1pc
<b>Accessory</b>	Operator's Manual:	1pc
	Charger:	1pc
	Rechargeable 300mAH Ni-MH battery (9V):	1pc

**Warranty** 18 months (6 months for probe)

**Dimension (L×W×H):** 142mm×112mm×42mm

**Color box Dimension (L×W×H):** 210mm×180mm×60mm

**Outer Carton Dimension (L×W×H):** 640mm×230mm×390mm

**Net Weight:** 300 g

**Gross Weight:** 470 g

**Gross Weight (20pcs):** 11 Kg