

# MEDIKRO primo

PC-based Office Spirometer

- ▶ Durable Office Spirometer
- ▶ More than 90 parameters
- ▶ Patient Safety
- ▶ User-Friendly



Product Code: M9492

Medikro® Primo is a professional and efficient PC-based spirometer for primary care. It is an accurate and easy-to-use system that complies with all the latest ATS/ERS spirometry standards.

Primo utilizes Medikro® Spirometry Software (MSS). The system provides general practitioners the assurance of reliable results with an intuitive navigation, adjustable interpretation options and clear visualization of results. With disposable Medikro® SpiroSafe flow transducers there is no need for bacterial filters or sterilization. This increases patient throughput and saves operational costs.



## Primo Highlights



### Patient safety

Disposable SpiroSafe flow transducers provide security in infection control. No exhaled air gets into the pressure tube or the device. This prevents patient-to-patient cross-contamination and protects the spirometer.


### User-friendly

Medikro Spirometry Software is user-friendly and provides intuitive navigation and use.



### Multiple calculation criteria and interpretation options

MSS provides multiple options for calculation of Best PEF, FEV% and FEF. In addition, interpretation options are available for ventilation function, spirometry findings and the bronchodilatation test.



### Trending

With the trending feature, doctors can track the historical development of one or more individual variables. This enables monitoring and documentation of treatment effects, identification of the need for treatment modifications, and evaluation of the changes in lung function.



### PDF export

Print final reports in PDF format or as images, or let the software do it for you automatically.

### Medikro eSupport

Medikro eSupport is your direct link to Medikro. It offers quick response time in case of support issues and provides online assistance 24/7. There is a vast knowledge base, an interactive troubleshooter, and an FAQ and download section. Live online support is available during office hours.



## Consumables & Accessories



**MEDIKRO® SpiroSafe**

Product code: M9256



**MEDIKRO® Calibration Syringe**

Product code: M9474 (3L) / M9477 (1L)



**MEDIKRO® Ambi**

Product code: M9489



🏠 Kellolahdentie 27, FI-70460 Kuopio, Finland  
✉ P.O. Box 54, FI-70101, Kuopio, Finland  
☎ +358 44 777 7002  
📞 +358 17 283 3300  
@ sales@medikro.com  
🌐 www.medikro.com

Supplied by

**team** medical supplies

✉ P.O. Box 6245 Dural DC, NSW 2158  
☎ 1300 22 44 50  
☎ 1300 22 44 60  
@ info@teammed.com.au  
🌐 www.teammed.com.au