


Cuticular


Reverse Cutting

Ethicon Syneture B. Braun

3/8 Circle



1/2 Circle



Straight Cutting




Plastic / Cosmetic

Premium Reverse Cutting with Micropoint

Ethicon Syneture B. Braun

3/8 Circle




Premium Cutting with Micropoint

3/8 Circle



1/2 Circle




Orthopaedic

Reverse Cutting

Ethicon Syneture B. Braun

1/2 Circle



1/2 Circle Heavy Bodied



Endoscopic


Taper Point

Ethicon Syneture B. Braun

Straight



1/2 Circle



B | BRAUN

SHARING EXPERTISE

Aesculap, Inc.

3773 Corporate Parkway
Center Valley, PA 18034
Phone: 800-282-9000


www.aesculapusa.com/sutures

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2007 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA. DOC571 REV. B 5.5M 7/07
Aesculap is an equal opportunity employer

B. Braun Closure Technologies

B. Braun Needle Comparison Chart



- Taper Point
- ▼ Reverse Cutting
- ▼ Lancet
- ⊕ Trocar Point
- Blunt Point
- ▲ Tapercut
- ⊙ Short Cutting Tip
- ▼ Reverse Cutting Micropoint
- ▼ Lancet Micropoint


B | BRAUN
SHARING EXPERTISE

General Surgery

Taper Point

Ethicon Syneture B. Braun


1/2 Circle



5/8 Circle



1/2 Circle Heavy Body



Blunt Point




MicroSurgery

Taper Point

Ethicon Syneture B. Braun

3/8 Circle



1/2 Circle




Cardiovascular


Taper Point

Ethicon Syneture B. Braun

3/8 Circle



1/2 Circle




Short Cutting

3/8 Circle



Trocar Point


1/2 Circle



Taper Cut

Ethicon Syneture B. Braun

1/2 Circle



3/8 Circle



Ophthalmic


Premium Spatula with Micropoint

Ethicon Syneture B. Braun

1/4 Curve



3/8 Circle



1/2 Circle




Progressive Curve



1/8 Circle



1/4 circle



Premium Reverse Cutting with Micropoint

3/8 Circle



1/2 Circle

