

## PRODUCT OVERVIEW

# Monosyn® Quick

Structure	Monofilament
Colour	Undyed
Chemical composition	Glyconate, a polymer made of: <ul style="list-style-type: none"><li>72 % Glycolide,</li><li>14 % Trimethylene carbonate</li><li>14 %-caprolactone</li></ul>
Sizes	USP 1 (metric 4) - 6/0 (metric 1)
Knot tensile strength retention	70 - 80 % after 5 days 50 % after 7 days 20 - 30 % after 10 days 0 % after 14 - 21 days
Mass absorption	Degradation by hydrolysis within 56 days
Sterilisation	Ethylene oxide

