


# Safil® Quick



Short-term synthetic absorbable braided and coated suture made of 100% polyglycolic acid (PGA)

### Advantages:

- High initial tensile strength
- Flexibility and outstanding handling properties
- Excellent knot security
- Easy gliding through tissue minimizes trauma



### Tensile Strength:

50% tensile strength at 7 days  
0% tensile strength at 14 days

### Absorption:

Mass absorption in approximately 42 days

### Surgical Indications:

Safil Quick is indicated for use in soft tissue approximation where short-term wound support is required and where rapid absorption would be beneficial.

### Common Applications:

- Subcutaneous and intracutaneous closure
- Gynecology
- Odontology
- Pediatric surgery
- Ophthalmic surgery
- Ligatures

