




Monosyn®



Monofilament synthetic mid-term absorbable suture made of glyconate

Advantages:

- High initial tensile strength
- Ideal degradation for soft tissue
- Extremely flexible
- Smooth passage through tissue



Tensile Strength:

50% tensile strength at 14 days
0% tensile strength at 28 days

Absorption:

Mass absorption in 60-90 days

Surgical Indications:

Monosyn is indicated for use in general soft tissue approximation and ligation where mid-term wound support is required.

Common Applications:

- Subcutaneous and intracutaneous closure
- Plastic and reconstructive surgery
- Gynecological surgery
- Gastrointestinal anastomosis
- Urology

