

# Gypsona® Extra 3

Fast Setting & Durable Plaster of Paris



# Gypsona®

# Gypsona® Extra 3 Plaster of Paris


Created for busy cast rooms and fracture clinics, Gypsona Extra 3 features a quick setting time for the fast immobilization of fractures. Once set, its composition ensures that the cast is strong and durable.

- Leno weave gauze for good conformability and stability
- Rich creamy texture
- Low plaster loss
- Moisture resistant packaging
- Rapid immobilisation of fractures

Gypsona® Extra 3 Plaster of Paris (POP) bandages are made from a leno weave gauze, uniformly impregnated with a high quality fast setting Plaster of Paris, which has a rich creamy texture for easy handling. The interlocked weave of the leno gauze makes the wet bandages stable, conformable and easy to mould as well as minimising plaster loss.

## Indications

- External splinting of fractures
- Immobilisation of diseased bones and joints
- Correction of deformities
- Emergency splinting
- Making negative and positive casts



## Ordering Codes

REF	Measurements	Units
65072-00	7.5cm x 3.5m	24 rolls/box
65072-01	10cm x 3.5m	24 rolls/box
65072-02	15cm x 3.5m	12 rolls/box
65072-03	20cm x 3.5m	12 rolls/box