

HYFRECATOR 2000 POWER SETTINGS



Procedure	Hyfrecator 2000
Telangiectasia	0-3 watts
Skin tag	1-3 watts
Mole or Nevi	2-6 watts
Seborrheic Keratoses	3-8 watts
Skin Cancer	2-10 watts
Venereal Warts	5-15 watts

ConMed Linvatec Australia Pty Ltd, 4/ 10 Rodborough Rd, Frenchs Forest NSW 2086

Supplied by

team medical
supplies

✉ P.O. Box 6245 Dural DC, NSW 2158

☎ 1300 22 44 50

📠 1300 22 44 60

@ info@teammed.com.au

🌐 www.teammed.com.au