



Transforming skin lesion treatment





Achieve exceptional treatment results

Freezen offers Precision, Practicality, and Effectiveness,
for a Rapid, No Stress treatment.

1

Forged from the widespread acceptance that cryotherapy is the benchmark treatment for all skin lesions.

2

Freezen provides a practical and easy to use instrument, that is cost effective and treats in seconds.

3

Specifically designed primarily for general practitioner, dermatology and cosmetic applications, it is the optimal solution for treating warts, condyloma, hemangiomas, skin tags, keratoses and all other benign skin disorders.

4

A no stress treatment. Freezen's high precision limits damage to healthy skin, minimizing discomfort to patients.

5

Effective result achieved from first treatment, with Freezen's high pressure achieving deeper penetration.





Practical Pen design - Optimal safety, no maintenance, no side effects. Uses the direct application of nitrous oxide at -89°C . It penetrates at 1mm per 5 seconds to a maximum depth of 5 mm.



Cost Effective Solution with low priced small cryogen cartridges.

« Our clients love it and so do we. We recouped the investment in just 30 treatments! »





Freezen is manufactured by Clinic6. We understand the need for opening new ways to improve the quality of healthcare delivery through our continuous dialogue with our medical markets. We monitor to change customer demands and market trends into practical solutions for everyday use.



Box with 15 N₂O capsules of 8g each.
Ref. CL-N2O 8G-BOX
Carton with 16 boxes.
Ref. CL-N2O-8G-CARTON



Freezen instrument set with user manual
and tweezers for filter placement.
Prepared for treatment of lesions Ø from 2-6mm.
Ref. CL-FP-00