HEM-907 Professional Blood Pressure Monitor











Pale Mounting Kit

Accurate and reliable

Complies with AAMI' standards for accuracy in automatic blood pressure measurement

Easy for operation and training

With advanced fuzzy logic technology, the HEM-907 determines ideal cuff inflation level based on patient arm size and blood pressure fluctuation. Simple push button operation and automatic deflation speed makes training clinical staff fast and easy

Dual measurement mode

AUTO Mode: the oscillometric measurement method

 Cuff automatically inflates to the optimal pressure according to patient's condition

MANU Mode: the auscultation measurement method, using a stethoscope

The unit measures blood pressure with the use of a stethoscope

Quiet, ideal for clinical environment

The HEM-907 operates quietly reducing stress on patients and clinical staff; perfect for night time use by your patient's bedside

Wall Mounting Kit

HEM 907 Accessories

•	Connector Plastic	DAVHEM90704	
	AC Adaptor	DAVHEM907ADAPT	
	Bladder Small	DAVHEM9078LDRS	
	Bladder Medium	DAVHEM9078LDRM	
	Bladder Large Adult	DAVHEM9078LDRL	
	Cuff Small 17 - 22cm	DAVHEM907CUFFS	
	Cuff Medium 22 - 32cm	DAVHEM907CUFFM	
	Cuff Large 32 - 42cm	DAVHEM907CUFFL	
	Air Tube 1.3 Mtr	DAVHEMTUBE130	
	Battery Pack	DAVHEM907PBAT	

Optional installation options				
Pole Mounting Kit	DAVHEM907PKIT			
Wall Mounting Kit	DAVHEM907WKIT			
Mobile Stand	DAVHEM907STAND			

Also available:

Extra Large Cuff & Bladder Set (XL) 42 - 50cm; includes 1m tubing DAVHEM907CUFFXL







Mobile Stand

Supplied by



№ P.O. Box 6245 Dural DC, NSW 2158

~ 1300 22 44 50

1300 22 44 60

@ info@teammed.com.au

mww.teammed.com.au