

# CryoClear®

Repair. Renew. Rejoice.

FROM BEFORE AND AFTER  
TO BEFORE AND **WOW.**




## The all new, clearly better cryotherapy solution

CryoClear® uses Carbon Dioxide which freezes at  $-79^{\circ}\text{C}$  to a depth of no more than 2mm into the skin, delivering a very safe and clinically effective treatment to superficial skin indications. CryoClear's high efficacy - over 95% - combined with less destruction of healthy tissue makes it perfect for treating age spots, sun spots, and skin tags. Whether used as a stand-alone treatment, or as a combination treatment to chemical peels, facials, micro-needling or dermaplaining, CryoClear offers a better, faster cryotherapy solution with significant benefits for clients and estheticians.

### SKIN TAGS | Treatment Time: One 5 second treatment (base of Stalk)



### AGE SPOT | Treatment Time: One 5 second treatment



Because CryoClear superficially treats indications that live within the epidermis, healthy skin will remodel to natural skin color as melanin is produced and lifted to the upper layers of the epidermis. This may take weeks to months following treatment.

[CryoClearRx.com](http://CryoClearRx.com)