



**Magnetically quick & simple**

DermLite DL4 can be magnetically attached to a DermLite iPhone 6 case without any additional parts.

- + compact all-aluminum design
- + now 45% brighter (DL4W only)
- + instant-on & auto shutoff
- + large 30 mm 10x lens system
- + charge from any USB port
- + polarized & non-polarized
- + PigmentBoost™ Plus (not in DL4W)
- + contact & non-contact dermoscopy
- + magnetically attach cameras, smartphones & iPad
- + IceCap™ infection control system
- + four charge indicator levels
- + two brightness levels
- + magnetic faceplate
- + sturdy protective belt holster

**IceCap™**

Disposable snap-on infection control



**The most refined vision in**

The 4th-generation DermLite dermoscope has been designed from the ground up to be more refined in every detail. An amazingly smooth, ergonomically optimized all-aluminum design that comfortably molds to your hand, DermLite DL4 puts phenomenal versatility, the largest optics of any pocket-size DermLite, and instant on/off control right at your fingertips. DL4 offers a larger field of view, and is still more than 25% smaller than the previous generation. DL4 uses magnetic connectivity to achieve groundbreaking convenience and fast workflows. Not only are both faceplate (with 10 mm reticle) and eyepiece magnetically attached, your smart devices (like



**dermoscopy.**

iPhone, iPad or Galaxy S series phones) also attach magnetically, enabling you to go between taking clinical and dermoscopic photos quicker than ever before. Along with the growing need for dermoscopy, infection control is becoming an increasingly important aspect of providing care. With its Infection Control and Epidemiology system IceCap™, DL4 has you covered. Disposable snap-on IceCaps attach to the spacer and may greatly reduce the risk of cross-contamination. A newly developed DL4 retractable glass faceplate design allows you to go between non-contact & contact modalities smoother than ever and lets you precisely focus your image while in skin contact.

DermLite DL4's high-powered LEDs and rechargeable battery produce light in an extremely energy-efficient way. And, in order to know when you run out of battery, DL4 now has a more granular four-level charge indicator. If you need to recharge, you can do so from any USB port using the standard Micro USB to USB cable included with each DL4.

**Versatile DL4 vs. Ultra-Bright DL4W**

While model no. DL4W produces ~45% more light from its exclusively white LEDs, the top-of-the-line model no. DL4 provides enhanced imaging of deeper pigmentation with two levels of color spectrum control (PigmentBoost Plus).