


# H.R.T TRAY for Hormone Replacement Therapy

Femcare Limited

**Ensuring a simple  
procedure remains  
simple...**



FEMCARE

# H.R.T Tray for Hormone Replacement Therapy

## Femcare Limited

The Femcare disposable sterile H.R.T Tray for hormone replacement therapy is designed to help the busy nurse and doctor reduce the time consuming preparation normally required prior to performing hormone implants. The H.R.T Tray eliminates both expensive re-sterilisation costs and also the time lost in cleaning instruments. The likelihood of delays caused by the occasional forgotten item is also eliminated, thus ensuring that a "simple procedure remains simple..."

The tray includes the necessary equipment for hormone implantation, which should be decided upon according to clinical indications.

### The Hormone Replacement Therapy Tray consists of two packs:-

The contents of Pack One are supplied sterile as listed below and are contained within a tray:-

#### Pack One

1 x Sterile plastic tray with reservoirs for antiseptic and hormone pellets  
1 x No. 11 Scalpel  
1 x Pair of Forceps  
1 x Pair of Surgeon Gloves  
5 x Cotton Wool Balls  
1 x Surgical Drape  
5 x Gauze Swabs  
2 x Hand Towels  
1 x Trocar (Obturator) and Cannula

Pack Two is an accessory to Pack One and contains individually packaged sterile items as listed below:-


#### Pack Two

1 x 5ml Syringe  
1 x 21G Hypodermic Draw-Up Needle (Green)  
1 x 23G Hypodermic Injection Needle (Blue)  
1 x Wound Dressing  
1 x Patient Record Card

#### Items Not Included

Local Anaesthetic  
Hormone Pellet

Product Code	Description
HT-WOS-853	HRT Tray, 1 each
HT-TCO-852	Trocar (Obturator) and Cannula, 10 each per box



## FEMCARE

### Europe

Femcare Ltd.  
32 Premier Way  
Romsey, Hampshire SO51 9DQ, UK  
Tel +44.1794.525.100  
Fax +44.1794.525.101  
enquiries@femcare.co.uk

### Canada

Femcare Canada  
6355 Kennedy Road, Unit 15  
Mississauga, ON L5T 2L5  
Tel +1 905.795.1102  
Fax +1 905.795.8154  
customer@femcare.ca

### Australia

Femcare-Australia Pty Ltd  
Unit 12, 5 Gladstone Road  
Castle Hill, NSW 2154  
Tel +61.2.9045.4100  
Fax +61.2.9894.1913  
customer@femcare.com.au