

**K** The "K" Series defines the HEINE-Standard full-size instrument for professionals, far beyond usual market standards. The line is characterised by unbreakable polycarbonate construction and full-feature performance to ensure the highest level of diagnostic accuracy.

### K180® F.O. Ootoscope

The standard, all around Fiber Optic Ootoscope

Ootoscope with Fiber Optic (F.O.) Illumination constructed of light, high-impact polycarbonate.

- Fiber Optic Illumination
- Polycarbonate construction
- XHL Xenon Halogen Technology
- Acrylic viewing window with 3 x magnification
- Swivelling viewing window built into instrument
- Insufflation port

### K180 Ophthalmoscope

A full-featured instrument at an economical price

- Precision HEINE optics. Brilliant retinal image
- Bright, white light. Clear retinal image
- Dust-resistant housing
- High-impact polycarbonate housing for durability
- 27 lenses from -35D to +40D
- Brilliant image even at high corrections

## HEINE K180® Diagnostic Set

The ultimate in performance and quality

Contents

- K180 Ophthalmoscope
- K180 Fiber Optic Ootoscope
- 1 set of reusable tips
- 1 hard case



HEA279101188

K180 diagnostic sets				
With 2.5V BETA battery handle			HEA279101188	*** \$996.57
With BETA4 USB rechargeable handle with USB cord and plug in power-supply (Li-ion battery)			HEA27927388	*** \$1400.57
With BETA4 NT rechargeable handle with NT4 charger (Li-ion battery)			HEA27923420	*** \$1792.01