

Blood-pressure cuffs and accessories

Latex-free blood-pressure cuffs

Cuffs for GAMMA G7/G5/GP/XXL LF sphygmomanometers.



- ✦ Latex-free.
- ✦ Fit the most widely used sphygmomanometers.
- ✦ High-quality velcro cuffs made of wear-resistant material. Physiologically tested.
- ✦ Double-sewn.
- ✦ Washable.
- ✦ Index mark for exact positioning.
- ✦ Clear indication of size and control mark for arm circumference.
- ✦ Extra thick, highly flexible tube. For a long life.
- ✦ Proven reliability. Absolutely airtight system.

	GAMMA G7/G5/GP	Single-tube cuffs (cpl.)			Single-tube bladders		Outer cuffs	
		Size	Circumference/cm		Size		Size	
	Infant [01]	6x28	8–13	M-000.09.700	5x15	M-000.09.661	6x28	M-000.09.059
	Child [02]	10x34	13–20	M-000.09.702	8x20	M-000.09.663	10x34	M-000.09.060
	Adult small [03]	14x53	20–29	M-000.09.704	11x23	M-000.09.664	14x53	M-000.09.061
	Adult [04]	14x58	29–41	M-000.09.619	12x28	M-000.09.665	14x58	M-000.09.055
	Adult large [05]	17x70	35–47	M-000.09.706	13x34	M-000.09.668	17x70	M-000.09.058
	Thigh [06]	20x86	>41	M-000.09.615	18x52	M-000.09.667	20x86	M-000.09.057

	GAMMA XXL LF	Twin-tube cuffs (cpl.)			Twin-tube bladders		Outer cuffs	
		Size	Circumference/cm		Size		Size	
	Infant [01]	6x28	8–13	M-000.09.701	5x15	M-000.09.671	6x28	M-000.09.059
	Child [02]	10x34	13–20	M-000.09.703	8x20	M-000.09.673	10x34	M-000.09.060
	Adult small [03]	14x53	20–29	M-000.09.705	11x23	M-000.09.674	14x53	M-000.09.061
	Adult [04]	14x58	29–41	M-000.09.629	12x28	M-000.09.675	14x58	M-000.09.055
	Adult large [05]	17x70	35–47	M-000.09.631	15x34	M-000.09.678	17x70	M-000.09.058
	Thigh [06]	20x86	>41	M-000.09.625	18x52	M-000.09.677	20x86	M-000.09.057

