



# Basics\* Tri-Layer AAMI2 Over-The-Head Isolation Gown



HALYARD\* AAMI Isolation Gowns provide barrier protection based on AAMI guidelines and helps reduce the risk of contamination and transmission of infectious organisms that lead to HAIs. HALYARD\* Isolation Gowns made from Tri-Layer fabric give clinicians the confidence and comfort they need to deliver the best possible patient care.

- Thumbhook with elastic cuff
- Closed back (with vent)
- Medium weight 3-layer SMS
- Meets AAMI Level 2 requirements
- Large size  
Gown Length: 127 cm  
Chest: 75 cm  
Length from Neck: 122 cm  
Sleeve length outer seam: 68 cm

**Code: HYH43008**

**Pkt/10**



**Looking after your team.** We are a proud independently owned Australian company

<b>Team Medical Supplies</b>	Phone 1300 22 44 50	info@teammed.com.au
	Fax 1300 22 44 60	www.teammed.com.au
<b>BP Medical (WA)</b>	Phone 1300 33 55 60	info@bpmedical.com.au
	Fax 1300 33 55 70	www.bpmedical.com.au

