



## Underpads - Full and Half Size

Absorbent & Versatile Protection for Patients and Surfaces

**Our  
Commitment  
to you**

**Better  
Infection  
Protection**

Looking for a product that's hygienic for patients and residents, as well as soft and absorbent for patient wellbeing? A product that will protect surfaces, linen and furniture is a variety of settings such as hospitals, aged care, General Practitioners, podiatry, veterinarian and obstetrics? A product that's versatile, waterproof and available in various absorbency levels?

### **Our Solution: Halyard\* Underpads - Full and Half Size**

HALYARD\* Underpads are constructed from a soft top layer and highly absorbent layers of tissue with a strong blue polyethylene backing that is quiet, drapable, and helps prevent strikethrough of liquid. Underpads are the ideal choice as a medical procedure underlay where fluids are involved, providing hygienic protection for patients and surfaces. HALYARD\* Underpads also assist in reducing the time taken for staff to clean and change linen.



## Application

- Underlay for wound care, dressing, procedures and surgical preparation
- Protection for linen and mattress
- Light incontinence aid
- Chair and surface protection
- Pillow protector
- Bedsheet
- Benchcover

## Dispensing

- Available in 5-ply: Full size and half-size
- Available in 8-ply: Full size

Inner poly bags- make it easier to store on shelves or on trolley, easier to divide up a shipper when smaller quantities are required.



## Absorbent tissue layers

Available in different sizes and in 5 ply/8ply to suit various absorbency needs.

## Strong polyethylene backing

Helps prevent strikethrough of liquid to keep surfaces dry.

## Soft top layer

Comfortable, gentle of patients' skin.

## Soft, smooth poly base

Quiet to assist with patient discretion and drapable to protect varied surfaces.

## Blue folded edges

Helps to contain fluid by preventing run off, blue edge is a visual guide for edge of pad against white bedsheet.

## Disposable

Hygienic, single use helps prevent cross contamination.

## Product Information

Code	Description	Size	Packaging
2765	HALYARD* Underpads 5 ply	39.5 x 55cm	4 x 75 sheets per polybag 300 sheets/case
49525-A	HALYARD* Half-Size Underpads 5 ply	39.5 x 27.5cm	6 x 100 sheets per polybag 600 sheets/case
2768	HALYARD* Underpads 8 ply	39.5 x 55cm	4 x 50 sheets per polybag 200 sheets/case



**O&M Halyard Australia Pty. Limited**

☎ 1800 664 227 ✉ CustomerService\_ANZ@owens-minor.com 🌐 [www.halyardhealth.com.au](http://www.halyardhealth.com.au)

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2025 O&M Halyard, Inc. All Rights Reserved.

O&M Halyard Australia Pty Ltd, Suite 1, Level 1, Trinita 3, 39 Delhi Road, North Ryde NSW 2113, Australia. <https://www.halyardhealth.com.au> COPY-05147