

MONOCRYL* (poliglecaprone 25) Suture

Product Characteristics	MONOCRYL* Suture
Physical Material	Poliglecaprone 25 (copolymer of glycolide and epsilon-caprolactone)
Construction	Monofilament
Absorbable / Non-absorbable	Absorbable
Coating (if applicable)	N/A
Colours	Undyed – beige and Dyed – violet
Available Size Range	2 through 6/0
Tensile Strength (Wound Support)	Medium Term (3-4 weeks)
Tensile Strength Retention in Tissue	Undyed - 50-60% of original strength remains at 1 week - 20-30% of original strength remains at 2 weeks Dyed - 60-70% of original strength remains at 1 week - 30-40% of original strength remains at 2 weeks
Absorption Time	Completely absorbed in 91-119 days

*Trademark