



STRONG & ABSORBENT WATERPROOF PROTECTION HALYARD* Protectapad

TRUSTED PROTECTION

FOR SURFACES & PATIENTS

SO YOU CAN CARE WITH CONFIDENCE



Looking for a product that helps to protect your patients' clothing as well as surfaces from all types of liquid such as blood, saliva, water, chemical seepage and spills? A product that's versatile, waterproof and absorbent?

OUR SOLUTION IS HALYARD* PROTECTAPAD

HALYARD* Protectapad is constructed from four highly absorbent layers of tissue and a layer of strong polyethylene backing, laminated together to help prevent strikethrough. The embossed pattern increases liquid absorbency, helping to prevent run off. Softness ensures optimum patient comfort and also assists with cushioning and protecting instruments. HALYARD* Protectapad is laminated for extra strength wet or dry and is single-use, making it the hygienic choice to protect patients and surfaces from staining and other spills.

APPLICATIONS

- Dental bib
- Hygienic surface for dental tools
- Medical procedures underlay
- Covering benches and trays
- Under IV arm of patients
- Pillow protector
- Surface liner for hairdresser/beautician products
- Patient bib
- Surface protector

DEPARTMENTS


- Dentists
- Oncology day units
- Physiotherapists
- GPs
- Optometrists
- Laboratories
- Beauty, hair & nail salons
- Pathology labs
- Aged care facilities
- Hospitals
- Blood banks

HALYARD* Protectapad

- ✓ **HIGHLY ABSORBENT TISSUE**
Absorbs fluid quickly
- ✓ **FOUR LAYERS OF SOFT TISSUE**
Provides thickness where you need it, for optimal patient comfort and to cushion and protect instruments
- ✓ **EMBOSSED PATTERN**
Increases liquid absorbency, helping to prevent run off
- ✓ **STRONG POLY BASE**
Provides good wet and dry strength¹ and helps to prevent strikethrough to keep patients and surfaces dry
- ✓ **DISPOSABLE**
Hygienic, helps reduce risk of cross-contamination

DISPENSING

HALYARD* Protectapad has a cut sheet presentation for easy, economical dispensing. And it's available in 2 sizes, small and large – so it helps better meet your needs. The small size contains 8 packs of 100 sheets, and the large size contains 4 packs of 100 sheets, with each individual pack wrapped in plastic for hygienic storage. This material is also available in a continuous roll: HALYARD* Benchroll (Code 2911).



PRODUCT INFORMATION

Code	Description	Size	Packaging
2706	HALYARD* Protectapad, Small	28.5cm x 21.5cm	100 sheets/pack 8 packs/case 800 sheets/case
2705	HALYARD* Protectapad, Large	28.5cm x 43cm	100 sheets/pack 4 packs/case 400 sheets/case

Our trusted range of products provide hygienic protection for surfaces and patients with a proven track record of consistent high quality, efficacy and performance.

SURFACES

HALYARD ISOWIPE* Bactericidal Wipes
 HALYARD VERSA* Towel
 HALYARD TERI* Wipers
 HALYARD TERI* Hydro Absorbent Wipers
 HALYARD* Absorbent Tray Liner
 HALYARD* Benchroll

PATIENTS

HALYARD* Underpads
 HALYARD* Bed Protector
 HALYARD* Protectapad
 HALYARD* Bedsheet
 HALYARD* Washcloth Extra
 HALYARD* Neutral Hand Wipes



Team Medical

Phone 1300 22 44 50 Fax 1300 22 44 60
 info@teammed.com.au
 www.teammed.com.au



BP Medical

Phone 1300 33 55 60 Fax 1300 33 55 70
 info@bpmedical.com.au
 www.bpmedical.com.au



O&M Halyard Australia Pty Ltd.

AUS: 1800 664 227

NZ: 0800 482 211

www.halyardhealth.com.au

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2018 O&M Halyard International. All Rights Reserved.

1. Mechanical break resistance test data on file. Note: These tests are directional based on small sample sizes.



Advantage*
HALYARD

Customers around the world choose Halyard Health because of our clinically-superior products, purposeful innovation, remarkable service and deep expertise.

It's what we call the HALYARD ADVANTAGE*