



# HIGHLY ABSORBENT & VERSATILE PROTECTION

## HALYARD\* Underpads

### TRUSTED PROTECTION

FOR SURFACES & PATIENTS

SO YOU CAN CARE WITH CONFIDENCE

Looking for a product that's hygienic for patients, as well as soft and absorbent for patient wellbeing? A product that will protect all types of surfaces, linen and furniture in a variety of settings such as hospital, aged care, GPs, podiatry, chiropractic, veterinarian and obstetrics? A product that's versatile, waterproof and available in various absorbency levels?

### OUR SOLUTION IS HALYARD\* UNDERPADS

HALYARD\* Underpads are constructed from a soft top layer and highly absorbent layers of tissue with a strong blue polyethylene backing that is quiet, drapable, and helps prevent strikethrough of liquid. Underpads are the ideal choice as a medical procedure underlay where fluids are involved, providing hygienic protection for patients and surfaces. HALYARD\* Underpads also assist in reducing the time taken for staff to clean and change linen.

### APPLICATIONS

- Protection for linen and mattress
- Bedsheet
- Light incontinence aid
- Chair and surface protection
- Pillow protector
- Underlay for dressings, procedures and surgical preparation
- Underlay for bedpan or cover to prevent splashing when removing
- Baby scale liner

### DEPARTMENTS

- Hospitals
- Aged Care Facilities
- Wards
- Obstetrics
- Outpatients
- Casualty
- Children's Wards
- Nursery
- Blood Banks
- GPs
- Allied Health
- Veterinary
- Early Childhood Clinics

# HALYARD\* Underpads

- ✓ **HIGHLY ABSORBENT<sup>1</sup> TISSUE LAYERS**  
Available in 5 ply and 8 ply to suit various absorbency needs
- ✓ **STRONG POLYETHYLENE BACKING**  
Helps prevent strikethrough of liquid to keep surfaces dry
- ✓ **SOFT TOP LAYER**  
Comfortable, gentle on patients' skin
- ✓ **SOFT, SMOOTH POLY BASE**  
Quiet to assist with patient discretion and drapable to protect varied surfaces
- ✓ **BLUE FOLDED EDGES**  
Helps to contain fluid by preventing run off, blue edge is a visual guide for edge of pad against white bedsheet
- ✓ **DISPOSABLE**  
Hygienic, helps prevent cross-contamination and also assists in reducing the time taken for staff to clean and change linen



## Team Medical

Phone 1300 22 44 50 Fax 1300 22 44 60  
info@teammed.com.au  
www.teammed.com.au



## BP Medical

Phone 1300 33 55 60 Fax 1300 33 55 70  
info@bpmedical.com.au  
www.bpmedical.com.au



## Advantage\* HALYARD

Customers around the world choose Halyard Health because of our clinically-superior products, purposeful innovation, remarkable service and deep expertise.

It's what we call the HALYARD ADVANTAGE\*

## DISPENSING

Available in two levels; 5 ply and 8 ply, HALYARD\* Underpads have a variety of applications including bed sheet / mattress protection, chair protection, pillow protection and bench cover.

- Inner poly bags - make it easier to store on shelves or on trolley, easier to divide up a shipper when smaller quantities are required. Each poly bag will have a label with a barcode, product code, product description and lot number allowing product to be identified easily outside the shipper. See below for packaging details.
- Easier to transport – no need to carry a large shipper, each poly bag can be transported throughout your facility.

## PRODUCT INFORMATION

Code	Description	Size	Packaging
2765	HALYARD* Underpads 5 ply	39.5cm x 55cm	4 x 75 sheets per polybag 300 sheets/case
2768	HALYARD* Underpads 8 ply	39.5cm x 55cm	4 x 50 sheets per polybag 200 sheets/case



**Our trusted range of products provide hygienic protection for surfaces and patients with a proven track record of consistent high quality, efficacy and performance.**

## SURFACES

HALYARD ISOWIPE\* Bactericidal Wipes  
HALYARD VERSA\* Towel  
HALYARD TERI\* Wipers  
HALYARD TERI\* Hydro Absorbent Wipers  
HALYARD\* Absorbent Tray Liner  
HALYARD\* Benchroll

## PATIENTS

HALYARD\* Underpads  
HALYARD\* Bed Protector  
HALYARD\* Protectapad  
HALYARD\* Bedsheet  
HALYARD\* Washcloth Extra  
HALYARD\* Neutral Hand Wipes



O&M Halyard Australia Pty Ltd.

AUS: 1800 664 227

[www.halyardhealth.com.au](http://www.halyardhealth.com.au)

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2019 O&M Halyard International. All Rights Reserved.

1. Absorbency test data on file.  
Note: These tests are directional based on small sample sizes.