



STRONG & ABSORBENT BENCH PROTECTION HALYARD* Benchroll

TRUSTED PROTECTION

FOR SURFACES & PATIENTS

SO YOU CAN CARE WITH CONFIDENCE



Looking for a product that helps to protect your surfaces from spills and staining, but is thick and soft enough to protect and cushion instruments and tools?

OUR SOLUTION IS HALYARD* BENCHROLL

HALYARD* Benchroll is constructed from four highly absorbent layers of tissue and one layer of strong polyethylene backing, laminated together to help prevent strikethrough. The embossed pattern increases liquid absorbency so spills are quickly soaked up, helping to prevent run off. Thick soft tissue ensures optimum cushioning, so it's suitable for instrument protection. It is laminated for extra strength wet or dry and is single-use, making it the ideal choice where fluids are involved, and the ideal way to protect bench surfaces from staining and other spills.

APPLICATIONS

- Bench protector
- Surface protector for dental or lab instruments
- Medical procedures underlay
- Pillow protector
- Linen protection

DEPARTMENTS

- Dentists
- Laboratories
- Pathology labs
- Hospitals
- Veterinary
- Allied health
- Fertility clinics

HALYARD* Benchroll

- ✓ **HIGHLY ABSORBENT TISSUE**
Absorbs fluid quickly
- ✓ **FOUR LAYERS OF SOFT TISSUE**
Provides thickness where you need it to cushion and protect instruments and equipment
- ✓ **EMBOSSED PATTERN**
Increases liquid absorbency, helping to prevent run off
- ✓ **STRONG POLY BASE**
Provides good wet and dry strength¹ and helps to prevent strikethrough to keep surfaces dry and hygienic
- ✓ **DISPOSABLE**
Hygienic, helps reduce risk of cross-contamination

DISPENSING

HALYARD* Benchroll is available as a continuous roll and can be used on any size bench by simply cutting the desired length, making it the economical choice for laboratories. This material is also available in a cut sheet format: HALYARD* Protectapad (Codes 2706 and 2705).

PRODUCT CONSTRUCTION & INFORMATION

| Code | Description | Size | Packaging |
|------|--------------------|--------------|---------------------------|
| 2911 | HALYARD* Benchroll | 41.5cm x 91m | 2 rolls/case 182m/case |

Our trusted range of products provide hygienic protection for surfaces and patients with a proven track record of consistent high quality, efficacy and performance.

SURFACES

HALYARD ISOWIPE* Bactericidal Wipes
 HALYARD VERSA* Towel
 HALYARD TERI* Wipers
 HALYARD TERI* Hydro Absorbent Wipers
 HALYARD* Absorbent Tray Liner
 HALYARD* Benchroll

PATIENTS

HALYARD* Underpads
 HALYARD* Bed Protector
 HALYARD* Protectapad
 HALYARD* Bedsheet
 HALYARD* Washcloth Extra
 HALYARD* Neutral Hand Wipes



Team Medical

Phone 1300 22 44 50 Fax 1300 22 44 60
 info@teammed.com.au
 www.teammed.com.au



BP Medical

Phone 1300 33 55 60 Fax 1300 33 55 70
 info@bpmedical.com.au
 www.bpmedical.com.au



Advantage* HALYARD

Customers around the world choose Halyard Health because of our clinically-superior products, purposeful innovation, remarkable service and deep expertise.

It's what we call the HALYARD ADVANTAGE*



O&M Halyard Australia Pty Ltd.

AUS: 1800 664 227

NZ: 0800 482 211

www.halyardhealth.com.au

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2018 O&M Halyard International. All Rights Reserved.

1. Mechanical break resistance test data on file. Note: These tests are directional based on small sample sizes