

# LIEBHERR

## Quality right down to the smallest detail

Liebherr vaccine and pharmacy refrigerators are Quality Care Pharmacy Program compliant and purpose-built with high quality materials to meet the stringent standards needed by your pharmacy. Trust in the superior performance and reliability of Liebherr to store your valuable vaccine and refrigerated medicines safely and effectively at the precise temperature-controlled conditions.



Precise temperature control. Vaccines and refrigerated medicines can be securely stored between +2 °C and +8 °C.





#### Integrated alarm system.

Audible and visual alarms alert users of temperature breaches, power failure and when door is left open for a set amount of time.

#### Integrated data memory.

Integrated data memory logs maximum and minimum internal temperature.



#### Temperature and alarm monitoring.

Alarms can be forwarded to remote monitoring systems using volt-free alarm contact. All vaccine and pharmacy refrigerators come with an external data logger for independent temperature monitoring and to comply with the "Strive for 5" Guidelines.



#### Battery back-up in event of power outage.

In case of power failure, the electronic controller of the LKPv 1420, LKPv 1423, LKPv 6520 and LKPv 6523 is powered by a 72 hour battery back-up which allows temperature and alarms to be monitored, and data to be transmitted to remote monitoring systems.



**Quality Care Pharmacy Program compliant.** Independently tested in accordance with the "Strive for 5" Guidelines and the QCPP Standard which requires temperatures to be stable between +2°C and +8°C.



### Maximum temperature consistency.

The forced-air cooling system and air flow provides a stable and consistent internal temperature, and highly efficient compressor offers good temperature recovery after door openings. Each unit is designed to EN 60068-3 for maximum temperature stability and reliable operation.



### Physical and keypad locks for security.

A physical lock protects against unauthorised access and keypad lock prevents temperature and alarm setting changes without a passcode.



### Versatile interior to meet dispensary demands.

The commercial grade polystyrol interior has channels to allow shelf height to be adjusted and optional Single AluCool<sup>®</sup> Drawers let users configure drawers to suit their individual needs.

### Vaccine & Pharmacy Refrigerators



### Superior Performance Refrigerators & Freezers

Gross capacity
Exterior dimensions in mm (w/d/h)
Interior dimensions in mm (w/d/h)
Energy consumption in 365 days <sup>1</sup>
Ambient temperature range
Refrigerant
Sound power level
Voltage / Connection rating
Cooling system
Defrost
Temperature range
Gradient <sup>2</sup> / max. fluctuation <sup>3</sup>
External cabinet finish
Door material
Internal liner material
Type of control
Temperature display
Power failure alarm
Fault: warning signal
Interface / Volt-free alarm contact
Interior light
Adjustable shelves

Interior light
Adjustable shelves
Usable shelf area in mm (w / d)
Shelf material
Shelf loading
Castors

Handle

Lock
Self-closing door
Door hinges
Gross / net weight

### Accessories

Shelf, plastic-coated
U-shaped trayslide, right
U-shaped trayslide, left
Foot pedal opener
Product temperature sensor
Roller Bars
Additional lock barrels (up to 10)

### LKPv 1423 MediLine

Medicine	Ivicui
1007 liter	1000 1
1397 litre	1366 li
1430 / 830 / 2160	1430/
1236 / 650 / 1500	1236 /
820 kWh	714 kV
+10 °C to +40 °C	+10 °C
R 290	R 290
58 dB(A)	58 dB(
220 – 240V~ / 4.0 A	220 -
forced-air	forced-
Automatic hot gas	Automa
0 °C to +16 °C	-2 °C
3.8 °C / 3.9 °C	3.6 °C
steel / white	steel /
glass door	steel
304-grade stainless steel	304-gi
electronic controller	electro
external digital	externa
immediately upon power failure for	immed
72 hrs	72 hrs
visual and audible	visual a
RS 485 / yes	RS 485
LED lighting, with separate switch	
8	8
1236 / 640	1236 /
wire shelves plastic-coated	wire sh
60 kg	60 kg
castors with brake at the front,	castors
castors at the rear	castors
full-length integrated handle	full-len
fitted	fitted
yes	yes
left hand / right hand	left har
267 / 224 kg	249/2
YL7112393	YL7112
YL9001761	YL900
YL9001757	YL900
1// 0500000	

### LKPv 1420 MediLine litre / 830 / 2160

1236 / 650 / 1550
714 kWh
+10 °C to +40 °C
R 290
58 dB(A)
220 – 240V~ / 4.0 A
forced-air
Automatic hot gas
-2 °C to +16 °C
3.6 °C / 2.9 °C
steel / white
steel
304-grade stainless steel
electronic controller
external digital
immediately upon power failure for
72 hrs
visual and audible
RS 485 / yes
8
1236 / 640
1236 / 640 wire shelves plastic-coated
wire shelves plastic-coated 60 kg
wire shelves plastic-coated 60 kg castors with brake at the front,
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear
wire shelves plastic-coated 60 kg castors with brake at the front,
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand 249 / 206 kg
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand 249 / 206 kg YL7112393
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand 249 / 206 kg YL7112393 YL9001761
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand 249 / 206 kg YL7112393 YL9001761 YL9001757
wire shelves plastic-coated 60 kg castors with brake at the front, castors at the rear full-length integrated handle fitted yes left hand / right hand 249 / 206 kg YL7112393 YL9001761

1 Measured at +25°C ambient temperature and +5°C set temperature for fridges and -20°C for freezers 2 Gradient as defined in EN 60068-3: the difference between the warmest and coldest average measurements, increased by their expanded uncertainty, throughout the duration. 3 Max. fluctuation as per EN 60068-3: greatest fluctuation value determined throughout the duration of measurements.

YL9590639 YL9590407

### Vaccine & Pharmacy Refrigerators



LKPv 6523	LKPv 6520	LKUv 1613	LKUv 1610
MediLine	MediLine	MediLine	MediLine
617 litre	602 litre	152 / 132 litre	142 / 133 litre
700 / 830 / 2160	700 / 830 / 2160	601 / 618 / 820	601 / 618 / 820
533 / 650 / 1500	533 / 650 / 1550	440 / 435 / 670	440 / 435 / 670
571 kWh	497 kWh	369 kWh	273 kWh
+10 °C to +40 °C	+10 °C to +40 °C	+10 °C to +35 °C	+10 °C to +35 °C
R 290	R 290	R 600a	R 600a
58 dB(A)	58 dB(A)	47 dB(A)	47 dB(A)
220 – 240V~ / 4.0 A	220 – 240V~ / 4.0 A	220 – 240V~ / 1.0 A	220 – 240V~ / 1.0 A
forced-air	forced-air	forced-air	forced-air
Automatic hot gas	Automatic hot gas	Automatic	Automatic
0 °C to +16 °C	-2 °C to +16 °C	+3 °C to +16 °C	+3 °C to +16 °C
2.0 °C / 1.6 °C	1.8 °C / 1.0 °C	5.1 °C / 4.9 °C	4.3 °C / 4.6 °C
steel / white	steel / white	steel / white	steel / white
glass door	steel	glass door	steel
304-grade stainless steel	304-grade stainless steel	commercial grade moulded polystyrol	commercial grade moulded polystyrol
electronic controller	electronic controller	electronic controller	electronic controller
external digital	external digital	external digital	external digital
immediately upon power failure	immediately upon power failure	when mains power return	when mains power return
for 72 hrs	for 72 hrs		i i i i i i i i i i i i i i i i i i i
visual and audible	visual and audible	visual and audible	visual and audible
RS 485 / yes	RS 485 / yes	RS 485 / yes	RS 485 / yes
LED lighting, with separate switch		LED lighting, with separate switch	
4	4	3	3
533 / 650	533 / 650	440 / 420	440 / 420
wire shelves plastic-coated	wire shelves plastic-coated	wire shelves plastic-coated	wire shelves plastic-coated
60 kg	60 kg	45 kg	45 kg
castors with brake at the front,	castors with brake at the front,		
fixed castors at the rear	fixed castors at the rear		
full-length integrated handle	full-length integrated handle	slimline handle with integrated	slimline handle with integrated
		opening mechanism	opening mechanism
fitted	fitted	fitted	fitted
yes	yes	yes	yes
right, reversible on site	right, reversible on site	right, reversible on site	right, reversible on site
164 / 138 kg	154 / 128 kg	45 / 43 kg	40 / 37 kg
YL7112393	YL7112393	YL7113333	YL7113333
YL9001761	YL9001761		
YL9001757	YL9001757		
YL9590659	YL9590659		
YL9590407	YL9590407	YL9590145	YL9590145
		YL9590521	YL9590521
		on request	on request

1 Measured at +25°C ambient temperature and +5°C set temperature for fridges and -20°C for freezers



Laboratory fridges	LKv 5710	LKv 3913	LKv 3910
with electronic controls	MediLine	MediLine	MediLine
Gross capacity	583 / 540 litre	386 / 332 litre	361 / 333 litre
Exterior dimensions in mm (w/d/h)	751 / 750 / 1844	601 / 618 / 1840	601 / 618 / 1840
Interior dimensions in mm (w/d/h)	634 / 538 / 1500	440 / 435 / 1635	440 / 435 / 1635
Energy consumption in 365 days <sup>1</sup>	438 kWh	480 kWh	309 kWh
Ambient temperature range	+10 °C to +35 °C	+10 °C to +35 °C	+10 °C to +35 °C
Refrigerant	R 600a	R 600a	R 600a
Sound power level	52 dB(A)	48 dB(A)	48 dB(A)
Voltage / Connection rating	220 – 240V~ / 1.5 A	220 – 240V~ / 1.5 A	220 – 240V~ / 1.5 A
Cooling system	forced-air	forced-air	forced-air
Defrost	automatic	automatic	automatic
Temperature range	$+3 \degree C to +16 \degree C$	$+3 \degree C$ to $+16 \degree C$	$+3 \degree C$ to $+16 \degree C$
Gradient <sup>2</sup> / max. fluctuation <sup>3</sup>	3.8 °C / 3.1 °C	5.6 °C / 4.7 °C	4.1 °C / 3.6 °C
External cabinet finish	steel / white	steel / white	steel / white
Door material	steel		steel
Internal liner material		glass door	
Type of control	commercial grade moulded polystyrol	commercial grade moulded polystyrol electronic controller	commercial grade moulded polystyr electronic controller
	electronic controller		
Temperature display	external digital	external digital	external digital
Power failure alarm	when mains power return	when mains power return	when mains power return
Fault: warning signal	visual and audible	visual and audible	visual and audible
Interface / Volt-free alarm contact	RS 485 / yes	RS 485 / yes	RS 485 / yes
Interior light	-	LED lighting, with separate switch	
Adjustable shelves	5	9	9
Usable shelf area in mm (w / d)	634 / 518	440 / 420	440 / 420
Shelf material	wire shelves plastic-coated	wire shelves plastic-coated	wire shelves plastic-coated
Shelf loading	60 kg	45 kg	45 kg
Castors	castors with brake at the front, fixed		
	castors at the rear		
Handle	slimline handle with integrated	slimline handle with integrated	slimline handle with integrated
	opening mechanism	opening mechanism	opening mechanism
Lock	fitted	fitted	fitted
Self-closing door	yes	yes	yes
Door hinges	right, reversible on site	right, reversible on site	right, reversible on site
Gross / net weight	99 / 89 kg	88 / 82 kg	70 / 65 kg
Accessories			
Shelf, plastic-coated	YL7113485	YL7113333	YL7113333
Adjustable legs		YL9590231	YL9590231
Additional lock barrels (up to 10)		on request	on request
Foot pedal opener	YL9094502	Unitequest	Uniequest
Roller bars	123034302	YL9590380	YL9590380
		15330300	12330300

1 Measured at +25°C ambient temperature and +5°C set temperature for fridges and -20°C for freezers 2 Gradient as defined in EN 60068-3: the difference between the warmest and coldest average measurements, increased by their expanded uncertainty, throughout the duration. 3 Max. fluctuation as per EN 60068-3: greatest fluctuation value determined throughout the duration of measurements.

### Shelves plastic-coated

Additional plastic-coated grid shelves can be retrofitted as required to provide extra storage options. The sturdy grid shelves have an extremely high loading capacity of up to 60 kg.

### Additional lock barrels

Up to 10 additional lock barrels to protect the stored products against unauthorised access are available as accessories for the laboratory appliances with electronic controller. In case of multiple appliances, different locks can be installed to ensure that the respective appliance can only be accessed by the staff members responsible for it.

### Foot-pedal door opener

Optional foot-pedal opener for hands-free use.

### Roller bars and adjustable legs

30 mm high roller sets can be retrofitted to the under-worktop LKUv models and height-adjustable legs can be retrofitted to the upright LKv models. The legs are adjustable between 115 mm and 170 mm in height.

In order to guarantee optimum stability, roller sets must not be fitted to appliances which are installed using a connecting frame. This is only permitted for individual appliances.







### Andi-Co Australia

1 Stamford Road, Oakleigh VIC 3166 Phone: 1800 685 899 Email: sales@andico.com.au www.andico.com.au

Whilst every effort is made to ensure all specifications and information in our brochure is correct at time of printing. Andi-Co Australia reserve the right to make changes without notification. Where actual appliance measurements, specifications and ratings are crucial we recommend that the actual appliance installation instructions that accompany each appliance be checked. Andi-Co Australia ABN 68 005 899 365 will not be liable for any loss arising from the information in this brochure, whether in negligence or otherwise. Images are for illustrative purposes only.