

Varicex® S/F

Zinc Paste Bandage

From A as in Application to Z as in Zinc

Varicex is a bandage saturated with zinc paste. It is ready for use and saves time when applying semi-rigid permanent bandages in thrombophlebitis, oedema treatment, chronic venous insufficiency, lower leg ulcers in the healing phase and follow-up care of fractures.

Varicex® S

Knitted bandage with two-way stretch for bandages applied using the circular technique; very easy to apply because of the stretch effect.

- viscose/polyamide
- two-way stretch (widthwise and lengthwise stretch)
- knitted bandage
- for bandages applied using the circular technique
- suitable for tapered limbs
- very easy to apply
- contains preservatives

Varicex® F

- cotton/viscose
- rigid bandage
- low moisture content for rapid dry time
- contains preservatives



Fields of application

- for applying semi-rigid permanent bandages in oedema treatment, chronic venous insufficiency, lower leg ulcers and thrombophlebitis
- for follow-up treatment of fractures

Properties

- bandages coated with zinc
- paste ready-for-use, time-saving
- non bulky
- selvedged edges
- white

Ordering information

Varicex® S

stretched length, individually packed in aluminium pouch

Width (cm x m)	REF	Shipping Units (pcs.)
10 x 5	20630	1/10/80
10 x 7	20631	1/10/60
10 x 10	20632	1/10/30

Varicex® F

unstretched length, individually packed in an aluminium pouch

Width (cm x m)	REF	Shipping Units (pcs.)
10 x 5	20600	1/10/80
10 x 7	20601	1/10/60
10 x 10	20602	1/10/50

For any queries please contact your local Lohmann & Rauscher (L&R) Sales Representative or L&R Customer Service on:

Australia: 1300 572 869

New Zealand: 0800 572 869

Varicex® is a registered trademark of Lohmann & Rauscher GmbH & Co. KG. © 2018

Lohmann & Rauscher Pty Ltd · 321 Ferntree Gully Road · Mount Waverley, Victoria 3149 · Australia
Lohmann & Rauscher Limited · 21 Queen Street · Auckland 1010 · New Zealand

www.Lohmann-Rauscher.com.au