








**LTR Medical**  
Building long-term relationships in healthcare

# Springfusor®

Spring-driven portable infusion pump



## Simple. Effective.

-  **Quality assurance**  
Ensures slow push is always given as a slow push
-  **Lightweight and portable**  
Allows patients to be ambulatory and requires no power source
-  **Flexible**  
Available in 10ml, 30ml or 50/60ml. Suitable for different flow rates.
-  **Cost-effective**  
Low acquisition cost
-  **Easy to use**  
No complex programming

# Springfusor<sup>®</sup> and Flow Control Tubing (FCT)

## Simple. Effective.

Engineered and calibrated using sodium chloride 0.9% to deliver a specific flow rate.



Springfusor<sup>®</sup> 10-28  
Springfusor<sup>®</sup> 30-AB  
Springfusor<sup>®</sup> 50/60

FCT are available in the following syringe volumes and flow rate times.\*

### SF10

- FCT 10ml in 5 minutes
- FCT 10ml in 15 minutes
- FCT 10ml in 30 minutes
- FCT 10ml in 60 minutes
- FCT 10ml in 4 hours
- FCT 10ml in 8 hours
- FCT 10ml in 24 hours

### SF30

- FCT 30ml in 5 minutes
- FCT 30ml in 15 minutes
- FCT 30ml in 30 minutes
- FCT 30ml in 60 minutes
- FCT 30ml in 2 hours
- FCT 30ml in 4 hours
- FCT 30ml in 24 hours
- FCT 30ml in 48 hours
- FCT 30ml in 96 hours

### SF50/60

- FCT 50/60ml in 15 minutes
- FCT 50/60ml in 30 minutes
- FCT 50/60ml in 60 minutes
- FCT 50/60ml in 2 hours
- FCT 50/60ml in 4 hours

\* Additional sizes and rates may be available on request.



W [ltrmedical.com](http://ltrmedical.com)

P 1800 319 419

E [info@ltrmedical.com](mailto:info@ltrmedical.com)



**Team Medical Supplies**  
P 1300 22 44 50 F 1300 22 44 60  
[info@teammed.com.au](mailto:info@teammed.com.au) [www.teammed.com.au](http://www.teammed.com.au)



**BP Medical (WA)**  
P 1300 33 55 60 F 1300 33 55 70  
[info@bpmedical.com.au](mailto:info@bpmedical.com.au) [www.bpmedical.com.au](http://www.bpmedical.com.au)