

UCS™ Debridement

A pre-moistened cleaning and debridement cloth



Ready to use debridement system to enable good cleansing of the wound, wound edge and care of the surrounding skin.

Suitable for

- Cleaning of acute, traumatic wounds
- Ulcer (edge and bed)
- 1st and 2nd degree burns
- Fistula abscesses
- Perilesional
- Interdigital
- General skin care

Product features

- Sterile and pre-moistened debridement cloth
- The cleansing solution contains poloxamer which helps to reduce the bio-burden, prevents skin infection and combats microbes¹
- Reduces bacterial colonization significantly²
- Removes hyperkeratosis and reduces scaling on the whole leg
- Non-allergic ingredients and a cloth made of viscose for a reduced risk of allergic reactions
- Keeps the wound area wet and elastic to support the healing process³
- No mechanical damage on the traumatic tissue³

Active ingredients	Effect
Aloe Vera	Makes the skin soft and supple; soothes and relaxes the skin
Allantoin	Increases the skin's moisture and helps to restore the natural defense function
Poloxamer	Helps to reduce bio-burden, prevents skin infection and combats microbes

Note

Go to:

www.medi.de/en/products/ucs-debridement to see an explanatory video, which demonstrates how to use the product.

Studies

- ¹ Downe A. How wound cleansing and debriding aids management and healing, *Journal of Community Nursing* 2014;28(4):33 -37.
- ² Cassino R. A New Dimension In Wound Cleansing, Poster Presentation 19th EWMA Conference, Helsinki 2009.
- ³ Palumbo FP et al. A new ulcer cleansing system in the management of „non-healing“ ulcers in elderly, Poster Presentation Tissue Viability Society, Kettering 2013.

Effective debridement, cleansing and moisturising system with high-quality ingredients

The solution within UCS Debridement will soften slough and necrotic tissue allowing the unique fibres of the device to catch and hold debris and effectively remove it from the wound bed. The skin friendly surfactant within the solution will disrupt and remove biofilm reducing the bacterial burden within the wound bed. The unique solution contains a keratolytic which will soften dry skin (hyperkeratosis) and in combination with the fibres of the cloth will loosen and remove these plaques.

Benefits for the clinician

- Unique, convenient and safe way to clean the leg and wound
- Ready to use
- One product for all wound types
- Both sides of the cloth can be used
- Can be combined with all modern wound dressings
- Treating the leg, treats the person
- Save time and improve skincare

Benefits for the patient

- Soft debridement cloth pre-moistened with a mild cleansing solution without damaging healthy cells
- Odour reduction
- Enables self-care
- Cleans and improves skin conditions with no need to soak
- Safe and simple to use



Article	Art.-No.
UCS Debridement (1 Box with 10 cloths) Size of one cloth: 20 x 20 cm	8972011



Manufacturer:
Welcare Research S.r.l.
Via San Giovanni sul Muro, 18
20121 Milano, Italy

Before / after comparison



Before UCS Debridement



5 minutes after using 1 x UCS Debridement