




ADVANCED
WOUND CARE

PLUROGEL®

BURN AND WOUND DRESSING



Micelle Matrix Technology™

A DISTINCT APPROACH TO WOUND MANAGEMENT



PluroGel is a burn and wound dressing utilising Micelle Matrix Technology. This is a concentrated, cell-friendly surfactant to aid in maintaining an optimal moist wound healing environment that promotes autolytic debridement. PluroGel helps protect healthy tissue and soften wound debris.

Promoting Patient Comfort

100% water soluble

PluroGel softens wound debris and allows it to easily rinse off during dressing change, creating a positive dressing change experience.

Gentle, non-irritating formula


This product facilitates a gentle experience for patients.

Better adherence to wounds

PluroGel maintains its consistency, enabling it to remain on the wound bed for better protection.



**MAINTAINS
CONSISTENCY**




PluroGel Micelle Matrix

PluroGel's micelle gel matrix maintains moisture in the wound and controls fluid loss, helping protect healthy tissue and soften wound debris.

Anatomy of a Micelle

ILLUSTRATION OF A MICELLE

PluroGel's micelle matrix contains surfactant micelles. A micelle's hydrophilic surface bonds with water, facilitating movement in the wound environment. Exudate and wound debris are absorbed within the hydrophobic core, allowing for easy removal during dressing change.



Ordering Information



Item No.	Item Description	Qty/Case
0102-170050	PluroGel 50 g jar	24/case
0102-170400	PluroGel 400 g jar	6/case

PluroGel PBWD (references 0102-170050 and 0102-170400) is a device in class IIa non-sterile.

These products are intended to be used by healthcare professionals.
Before use, consult instructions and precautions on the corresponding labelling.

**Contact your Medline Account Manager for further information
1800 110 511**



Medline Australia & New Zealand
2 Fairview Place, Marsden Park
Sydney NSW 2765
1800 110 511
www.medline.com/au

Supplied by:

team medical supplies

