

JX820D-1 EMERGENCY SUCTION PUMP

Emergency Portable Suction Series



MAINS OR BATTERY OPERATED SUCTION PUMP.



OFFICE AND SHOWROOM:

23/9-15 Sinclair St, Arundel QLD 4214

Phone: (07) 5529 0753 or Free Call: 1800 677712

General enquiries: sales@medeleq.com.au

Shop online now www.medeleq.com.au/Store

Distribution all States

Portable Suction Unit **JX820D-1**

- **Modern Designed A.B.S outer shell for easy maintenance. Simple controls for easy operation.**
- **A great balance between size, style and performance. This unit is light weight and produces low noise, but it still offers great flow rates and maximum suction pressure to ensure strong reliable suction required by Hospital, Clinic and other Medical Users.**
- **This unit incorporates a Twin Diaphragm pump, using aluminium alloy die cast body and various precision elements. High running hardness, wear resistance, non-polluting, long life and free maintenance. No lubrication or daily maintenance is required.**

Specifications

1. **Max Vacuum: =0.08MPa (600mmHg)**
2. **Vacuum range : 0.02MPa ~ Max vacuum**
3. **Flow rate : =20 L/min**
4. **Noise :=65dB(A)**
5. **Jar capacity: 1000mL**
6. **Operating voltage: DC 12V or AC 100V~240V 50/60Hz**
8. **JX820DBattery duration on full charge: =30min;**
JX820D-1Battery duration on full charge: =60min
9. **JX820D/JX820D-1Dimension: 410mm×205mm×420mm**
- 10.**JX820D:Net Weight 4.5kg**



Optional Carry Bag



Twin Piston,
Diaphragm Pump



Disposable Liners or
Reusable Collection Jar