INSTILLAGEL®

STERILE LUBRICANT lidocaine 2% available with and without chlorhexidine



INSTILLAGEL® LIDO

- Pain relief through the anaesthetising effect of 2% Lidocaine Hydrochloride.
- Excellent lubrication properties of hydroxymethylcellulose base adheres especially well to the urethra.
- Long-lasting lubrication due to the nondrying effects of added propylene glycol.
- 11ml syringe for that extra bit of lubrication.
- Now available as a certified and registered medical device.

Instillagel® Almed Sterile Lubricant.

INSTILLAGEL® ALMED

- All the benefits of Instillagel® Lido PLUS 0.05% Chlorhexidine for patients without Chlorhexidine sensitivity.
- Two added hydroxybenzoates clinically analysed to have extra efficacy against bacteria.
- Superior adhesion to the mucosa provides an antiseptic shield before, during and after the procedure.
- Clinically evaluated effective against MRSA.
- Closely monitored low allergy history over 48 years.

INSTILLAGE STERILE LUBRICANT lidocaine 2% available with and without chlorhexidine



BETTER FORMULA TO SUIT YOUR USE

Whether it's the added hydroxybenzoates clinically analysed to have extra efficacy against bacteria with the Chlorhexidine Almed version or the added propylene glycol to maintain non-drying effects in our standard Lido syringe hydroxymethylcellulose base, we have the sterile anaesthetic lubricant suited to your patient.

In the simple test of thickness, our formula is less viscous; designed to spread more easily up the urethra and coat the entire surface.

It is designed to adhere especially well to the mucosa and last through a standard transurethral resection of over 2000 resectoscope movements. It contains no fragment particles or micro bubbles and makes an excellent ultrasound medium.



STERILE PACKAGING THAT CONFORMS TO ASEPTIC TECHNIQUE

Instillagel® comes in an autoclaved blister pack for aseptic technique procedures without any assembly of the syringe prior to use. This conforms to ACORN guidelines for maintaining a sterile field:

" All items introduced onto a sterile field should be opened, dispensed and transferred by methods that maintain sterility and integrity. Edges of all sterile wrappers and packages are considered nonsterile. Therefore, when opening an item, be sure all the edges are secured so that they do not contaminate the field or the sterile packaged item."



SYRINGE TIP SPECIFICALLY DESIGNED FOR LESS TRAUMA ON ENTRY

The Instillagel® tip is designed to fit easily into the urethra with the base of the syringe acting as a stopper against backflow during delivery. The Instillagel® tip is highly polished and smooth unlike other product tips which can feel sharp and rough due to being produced using a mould. It won't cause damage to the meatus or cause unnecessary discomfort when the syringe is inserted and positioned correctly.

For patient safety, it's important to know that the tip of the Instillagel® syringe has also been designed not to fit a needle so there is no risk of injecting the solution by mistake.



Distributed in Australia by: InterPharma Pty Ltd Suite 103 / 39 East Esplanade, Manly NSW 2095 Ph: 02 9976 6876 www.interpharma.com.au Distributed in NZ by: Healthcare Logistics

58 Richard Pearse Drive, Airport Oaks, Mangere 2022 Ph: 02 9976 6876