

HCFM2 - Type IIR Disposable 3-Ply Face Mask

- High Bacterial Filtration Efficiency (BFE) $\geq 98\%$
- Protection against fluid exposure at a pressure level of 120 mmHg
- 100% Hypoallergenic
- Size: 175 x 95mm
- Colour: light blue
- Latex free
- 50 pieces per package
- Covers nose, mouth and chin
- To be used by general public and medical staff.
Not suitable for surgical environment.
- This 3-ply protects against infectious pathogens, pollen and dust particles -outer layer (blue) filters large particles, middle layer filters pathogens and the inner layer (white) provides comfort to the skin



50
PCS

