


## Laryngeal Mask Airways

A smooth one-piece design allows for easy insertion and guide markings assist with positioning. Any condensation, blood secretions and regurgitated gastric contents are easily visualised through the translucent tubing.

- ⊗ Available in adult, paediatric and neonate versions
- ⊗ Colour coded for simple size selection
- ⊗ Clear guide markers for easy positioning
- ⊗ Soft seal cuff reduces pressure and minimises potential trauma
- ⊗ A pilot balloon allows for visual and tactile assessment of the inflated / deflated cuff
- ⊗ 100% medical grade silicone
- ⊗ Translucent tubing for clear visual monitoring



Laryngeal Mask Airways				STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY
AN030000	Silicone Disposable Laryngeal Airway Mask	1	Purple	10
AN030001	Silicone Disposable Laryngeal Airway Mask	1.5	White	10
AN030002	Silicone Disposable Laryngeal Airway Mask	2	Red	10
AN030003	Silicone Disposable Laryngeal Airway Mask	2.5	Orange	10
AN030004	Silicone Disposable Laryngeal Airway Mask	3	Green	10
AN030005	Silicone Disposable Laryngeal Airway Mask	4	Yellow	10
AN030006	Silicone Disposable Laryngeal Airway Mask	5	Pink	10