

Infusion Set with NeutralSite™

DESCRIPTION

The Infusion Set includes a vented spike attached to a drip chamber with 15µm filter, roller regulator and 1 x NeutralSite™.

FEATURES & BENEFITS

- Vented perforator
- Male luer lock with rotating collar
- NeutralSite™
- Rigid roller regulator for accurate flow
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter 3.0mm
- Tubing length 220cm
- Priming volume: 18.0mL
- Drip chamber: 20 drops /mL
- Filter: 15µm
- NeutralSite™: Neutral displacement needleless access site
- Sterile EO
- ARTG Number: 304860
- The product is MRI and CT suitable
- Lipid Resistant: YES
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Not made with Natural Rubber Latex
- Material: LDPE, ABS, PVC, Nylon filter
- NeutralSite™ Materials: PP, ABS and Silicone
- Shelf Life: 3 Years

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- EN ISO 594-1 Conical fittings with a 6% [Luer] taper for syringes, needles and certain other medical equipment —Part 1: General requirements
- EN ISO 594-2 Conical fittings with 6 % [Luer] taper for syringes, needles and certain other medical equipment —Part 2:Lock fittings

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Drip Chamber
2	Tubing
3	Roller Clamp
4	NeutralSite™
5	Tubing
6	Rotating Male Luer
7	Cap

Infusion Set			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV010003	Infusion Set	220cm	50