

Urine Bag T-Tap

DESCRIPTION

The Urine Bags are used to collect and contain urine via an indwelling catheter or leg bag. Our range of Urine Bags are available with or without a non-return valve. Variations in tube length ensure the bag is optimised for different settings.

FEATURES & BENEFITS

- Durable material
- Clear volume graduation
- With or without non-return valve
- With T-Tap
- With or without needleless sampling port
- Varying tube lengths to support patient's needs


PRODUCT DATA

- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Material: PVC, PP, PE
- Shelf life: 5 years
- ARTG number: 267875
- Sterile EO

ISO STANDARDS

- ISO 14971:2012 Medical Device - Applications of Risk Management.
- ISO 10993-1:2009 Biological compatibility test report
- ISO 8669:1996 Urine collection bags

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	T-Tap
2	Bag
3	Tube
4	Connector
5	Cap

Urine Bag T-Tap

Urine Bags				STERILE
RE-ORDER	DESCRIPTION	TUBING LENGTH	BAG CAPACITY	CARTON/QTY
UR062000	Urine Bag T-Tap NRV ID 5mm OD 6.3mm	90cm	2000mL	250pcs
UR062001	Urine Bag T-Tap ID 5mm OD 6.3mm	90cm	2000mL	250pcs
UR062003	Urine Bag T-Tap NRV, Sample Port, Clear Tube, ID 7.4mm OD 9.8mm	120cm	2000mL	100pcs
UR062007	Urine bag 2000mL T-TAP NRV	110cm	2000mL	250pcs
UR062008	Urine Bag T-Tap	110cm	2000mL	250pcs