

Leg Bag T-Tap with Silicone Lined Straps

PRODUCT DATA SHEET



DESCRIPTION

The Leg Bags are designed to provide discreet storage of urine under clothing. They are attached to the leg with silicone lined straps to prevent the bag from sliding down the leg. The bags are made from non-woven material, soft to the skin. Available in; 350mL, 500mL and 750mL versions with 10cm, 30cm or 50cm tube with bonded stepped connector, anti-reflux and T-Tap.

FEATURES AND BENEFITS

- Discreet underneath clothing
- Easy to use T-Tap to empty bag
- Available in 350mL, 500mL & 750mL bag sizes
- Bonded stepped connector to tubing
- Anti-reflux valve
- Tubing lengths: 10cm, 30cm & 50cm
- Comfortable silicone lined straps for non-slipping
- Soft non-woven backing

PRODUCT DATA

- Latex free
- Material: PVC, PP, PE & Silicone
- Packaging: Paper and PE
- Straps: Urethane elastic fibre
- DEHP: Yes
- Sterile: EO
- Shelf life: 5 years
- ARTG Number: 267876

RE-ORDER	DESCRIPTION	TUBING LENGTH	CAPACITY	CARTON/QTY
UR135010	LEG BAG - T-Tap NRV Short Tube with Silicone Lined Straps	10cm	350mL	80
UR135030	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	30cm	350mL	80
UR135050	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	50cm	350mL	80
UR150010	LEG BAG - T-Tap NRV Short Tube with Silicone Lined Straps	10cm	500mL	80
UR150030	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	30cm	500mL	80
UR150050	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	50cm	500mL	80
UR175010	LEG BAG - T-Tap NRV Short Tube with Silicone lined Straps	10cm	750mL	80
UR175030	LEG BAG - T-Tap NRV Long Tube with Silicone lined Straps	30cm	750mL	80
UR175050	LEG BAG - T-Tap NRV Long Tube with Silicone lined Straps	50cm	750mL	80